

SIEMENS



www.bestsound-technology.com



Siemens Hearing Instruments. Product Portfolio – Spring/Summer 2015



Life sounds brilliant.

Siemens Hearing Instruments Product Overview.

We offer a broad range of hearing aids in many different models and performance levels to help Hearing Care Professionals meet their clients' needs, preferences and demands. Take advantage of this diversity to enhance your customers' satisfaction.

| Model | Performance levels | RIC | BTE LifeTube/ ThinTube | BTE Earhook | ITE IIC, CIC, ITC, ITE | Connexx |
|------------------------------|-----------------------------|---|---|--|---|---------|
| Ace binax | 7bx 5bx ^{New!} 3bx |  | | | | 7 |
| Pure binax | 7bx 5bx ^{New!} 3bx |  | | | | 7 |
| Carat binax | 7bx 5bx ^{New!} 3bx |  | | | | 7 |
| ^{New!} Motion binax | 7bx 5bx 3bx | |  |  | | 7 |
| ^{New!} Insio binax | 7bx 5bx 3bx | | | |  | 7 |
| Pure micon | 7mi 5mi 3mi |  | | | | 7 |
| Siemens Life micon | 7mi 5mi 3mi | |  |  | | 7 |
| Motion micon | 7mi 5mi 3mi | |  |  | | 7 |
| Aquaris micon | 7mi 5mi 3mi | |  |  | | 7 |
| Nitro micon | 7mi 3mi | | |  |  | 7 |
| Insio micon | 7mi 5mi 3mi | | | |  | 7 |
| Orion | Essential |  |  |  |  | 7 |
| Sirion | Essential | |  |  |  | 7 |
| Intuis | Basic | |  |  |  | 6 |
| Lotus | Basic | |  |  |  | 6 |

Always ahead.

Siemens hearing instruments have shaped and established a long history of innovation. For over 130 years now, new technological benchmarks have constantly been set. And this success story continues today: With binax™, the next generation of BestSound™ Technology, we have found a way to effectively support binaural hearing. Plus, the powerful binax platform is the world's first technology clinically proven to enable wearers to hear better than people with normal hearing in certain demanding environments.*

Expertise meets innovation.

With Insio binax™, we have taken this groundbreaking hearing aid technology one tremendous step further. Thanks to the revolutionary binaural OneMic directionality, now for the first time ever, hearing instruments with only one microphone can offer enhanced natural directional hearing. As a result, even wearers of tiny CICs can enjoy binaural directionality for a clearer speech understanding to the front even in demanding environments. A true world's first. And to ensure that our binax portfolio covers nearly all hearing aid wearers' needs and wishes, we have extended it even further by introducing Motion binax™. Motion SX and SA binax come in a sleek new design that hides completely behind wearers' ears. These BTE models offer supreme directionality, the utmost discretion and flexibility when used together with the new ThinTube. And the new Motion PX binax offers wearers a high level of amplification, wearing comfort, and handling convenience.

Added value for all.

Let us help you provide the best options for your customers and highlight your expertise. The extensive Siemens portfolio of products and professional hearing care support can give you the edge you need to succeed. Both you and your customers will benefit from the clarity of Siemens hearing instruments and the convenience of our technology.

We give our very best so we can achieve more together!


Contents

| | |
|---|----|
| History | 4 |
| binax BestSound Technology | 6 |
| What's new? | 10 |
| Hearing Solutions for Everyone | 12 |
| binax: Feature Overview | 14 |
| binax: Advantages of a Higher Performance Level | 15 |
| micon: Feature Overview | 16 |
| micon: Advantages of a Higher Performance Level | 17 |
| Ace binax | 18 |
| Pure binax | 20 |
| Carat binax | 22 |
| Motion binax | 24 |
| Insio binax | 26 |
| Pure micon | 28 |
| Siemens Life micon | 30 |
| Motion micon | 32 |
| Aquaris micon | 34 |
| Nitro micon | 36 |
| Insio micon | 38 |
| Orion | 40 |
| Sirion | 42 |
| Intuis | 44 |
| Lotus | 48 |
| SecureEar | 52 |
| Direct audio input and child-lock battery door | 54 |
| Accessories | 56 |
| easyTek and easyTek App | 58 |
| Easy pairing | 59 |
| touchControl App | 60 |
| App download | 61 |
| Audio Streamers | 62 |
| Remote Controls and eCharger | 63 |
| miniReceiver 2.0 | 64 |
| miniReceiver Fitting | 65 |
| Click Domes and Click Molds | 66 |
| ThinTube | 68 |
| Earmolds | 71 |
| Counseling Suite | 72 |
| Connexx 7 | 74 |
| Unity 3 | 75 |
| Workflow and Solutions | 76 |
| eARena | 77 |
| Point of Sale Support | 78 |
| Jewel Cases | 81 |
| All Sizes at a Glance – RIC and BTE | 82 |
| Detailed Feature Overview RIC | 84 |
| Detailed Feature Overview BTE | 86 |
| Detailed Feature and Color Overview Custom | 90 |
| Color Overview RIC and BTE | 94 |

* See legal note on final page in this document.

Siemens hearing aids have been helping the world to hear – for more than 130 years.

We have shown a particularly high level of commitment in the field of audiology. With one major goal: We would like to contribute our solutions to better hearing and increase the quality of life for people with hearing impairments. That was the case back in the times of Werner von Siemens, and it is still the case today.



“Personal achievements can only be recognized by the benefits they offer others.”

Werner von Siemens,
inventor of the first hearing instrument

1878

Werner von Siemens develops a telephone receiver that also helps the hearing impaired.

1913

Esha-Phonophor, Siemens begins the serial production of hearing aids.

1949

Phonophor Alpha, Siemens presents their pocket-sized hearing instrument.

1959

Auricullete 326, Siemens launches their first behind-the-ear hearing aid.

1966

Siretta 339, Siemens presents the first in-the-ear hearing instrument.

1987

Telos, Siemens launches the first remote control in the world of hearing technology.

1997

Prisma, Siemens presents the first digital hearing aid featuring two microphones.

2002

Triano, Siemens launches the first hearing instrument with three microphones.

2004

Acuris, Siemens introduces the first hearing instrument featuring a wireless system, e2e wireless, small enough to fit in CICs.

2005

iScan, Siemens introduces the first in-office digital ear impression scanner.

2008

Tek, Siemens introduces a revolutionary wireless enhancement system giving wearers unprecedented access to all favorite audio devices.

2010

BestSound Technology, Siemens writes a new chapter in hearing aid technology.

Life sounds brilliant.

2011

Aquaris, Siemens introduces the first digital waterproof, dustproof and shock-resistant hearing instrument.

2012

Siemens researchers win the “**Deutscher Zukunftspreis**” (German Future Award).

2014

With **binax**, the next generation of BestSound Technology, Siemens introduces a binaural hearing aid technology better than normal hearing.*

* See legal note on final page in this document.

2015

A world's first: The pioneering **binaural OneMic directionality** enables binaural directionality even in hearing aids with one microphone such as CICs.



World's First!

Now offering binaural OneMic directionality.

A world's first dedicated to binax custom hearing aids.

With binax, the world's first technology that outperforms normal hearing*, we raised the standard of Soundability™, drawing inspiration from the healthy auditory system. This next generation of BestSound™ Technology simulates the effects of binaural hearing and brings its advantages back to the wearer. It is so powerful that for the first time, wearers with a pair of binax models Carat, Carat A or Pure can understand speech in demanding environments even better than people with normal hearing.* The new wireless, two-microphone Motion binax and Inzio binax models have the same directional performance as Pure binax™ and Carat binax™.

* See legal note on final page in this document.

› Binaural OneMic directionality for hearing aids with one microphone: Yet another world's first from binax.

Until now, automatic directionality could not be offered for hearing aids equipped with only one microphone. But thanks to binax, now even very small custom hearing aids – such as Inzio binax CICs – are able to offer binaural OneMic directionality because in a bilateral setup, they can use the microphone signals of both instruments. As a result, even CIC wearers can enjoy the benefits of automatic directionality in noisy situations. And just like all other binax key features, binaural OneMic directionality works fully automatically with minimal power consumption.

› New features for all binax products.

binax continues to make unprecedented contributions to audiology. The outstanding new features "Adaptive streaming volume" and "Ocean wave tinnitus therapy signals" complete the already existing feature set of binaxFocus, binaxSound, and binaxGuide. And the best part is, they are also available for our existing binax models.




High Definition Sound Resolution for the best binaural listening experience ever.

In order to emulate natural binaural listening, we first need to link the two hearing aids, just like the brain uses input from both ears. Building upon the original e2e technology, **e2e wireless™ 3.0** enables exchange of audio signals between two hearing aids in a bilateral fitting. As a result, each hearing aid not only uses the information picked up by itself, but also the acoustic signals received by the contralateral hearing aid. In two-microphone instruments, this means each instrument works with input from four microphones. In a bilateral fitting, they create a virtual eight-microphone network. Similarly, e2e wireless 3.0 allows single microphone instruments such as the tiny Inzio binax CICs to exchange audio signals between each other. In bilateral fittings, they create a virtual four-microphone network. These virtual microphone networks provide a heightened sensitivity of the acoustic environment: We call this "High Definition Sound Resolution" or HDSR.


› binax facts

- e2e wireless 3.0 transmits 1,000 x more data per second than the previous generation
- 16 % more transistors than micon
- 20 % more instructions per second than micon
- 48 channels, 12 kHz bandwidth
- Binaural audio transmission technology with the industry's lowest power consumption

e2e wireless 3.0 virtually doubles the microphones and creates HDSR.



For binax hearing aids with two microphones, e2e wireless 3.0 creates a **virtual eight-microphone network**.



For binax hearing aids with one microphone, e2e wireless 3.0 creates a **virtual four-microphone network**.

binaxSound binaxGuide binaxFocus



› Advanced binax features for a new level of Soundability.


binaxFocus, binaxSound, and binaxGuide create the perfect equilibrium between sound quality and audibility for each individual wearer – raising the bar for Soundability.

Find out more about the exceptional binax features on the following pages.

binaxFocus.


Hearing that's more precise than ever – in all situations.

› Binaural OneMic directionality




As another world's first, binaural OneMic directionality lets even wearers of single-microphone hearing aids experience the benefits of automatic directionality – so they can follow conversations in noisy environments more easily.

› Narrow directionality



binaxFocus lets wearers focus precisely on their desired sound source, whether it comes from the front, back or side. To the front, it can narrow the beam of the directional microphones, so wearers can understand speech in difficult listening situations even better than people with normal hearing.*

› Spatial SpeechFocus




* See legal note on final page in this document.

binaxSound.


Features in harmony for the best audio output.

› Adaptive streaming volume New!



Adaptive streaming volume enhances the listening experience when streaming audio with easyTek™. No matter which situation, adaptive streaming volume automatically adapts the volume of the streamed audio according to the changing background noise level. This way, streamed audio is always audible above the noise. The streaming volume no longer needs to be manually readjusted, leaving wearers free to simply enjoy what they want to hear.

› eWindScreen binaural




binaxSound ensures optimal sound quality, even in challenging situations. In the presence of wind, eWindScreen™ binaural strategically transmits audio signals from the hearing aid with the better audio quality to the other ear. This effectively reduces the annoying effects of wind without compromising spatial perception or speech intelligibility.


binaxGuide.

Hearing has never been so individual.

› Spatial Configurator: Span




› Spatial Configurator: Direction



As an addition to automatic functions that adapt to changing acoustic environments, Spatial Configurator empowers wearers to control their hearing aids in unique listening situations. Via the easyTek and easyTek App, they can individually adjust the span and direction of their directional microphone focus beams, and decide how much of their environment they want to hear.

› Expanded tinnitus function with ocean wave therapy signals New!



The expanded tinnitus function in binax hearing aids is now more comprehensive than ever. In addition to the proven static therapy signals that offer preset or customizable solutions, wearers can now also choose from four nature-inspired ocean wave therapy signals. These therapy signals provide a positive, soothing, and stress-relieving listening experience ideal for tinnitus patients who prefer a more natural therapy signal to the artificial static ones.

Minimized design, maximized performance: Insio binax customized hearing aids.


World's First!

With the groundbreaking binax technology, we have set another benchmark in audiology. Until today, wearers of traditional custom hearing aids always had to decide between technology and size: Tiny, discreet hearing aids could not be equipped with many technological features, while high-tech custom instruments used to sit visibly in the ear. Now with Insio binax customized hearing aids, we managed to narrow this gap.


Smart construction and audio data exchange even in CICs.

By rotating the receiver, we gained space to optimally position the interior components. As a result, we could move an improved wireless coil into that space for a more compact design. An electromagnetic shield protects the components against interference. Thanks to the advanced wireless coil, now even CICs can be equipped with e2e wireless 3.0 and, for the first time ever, exchange audio data. As a result, custom hearing aids now can be built even smaller than ever before, while, at the same time, offering enhanced functionality.


New design with advanced coil



Hearing aids with one microphone each



Binaural OneMic directionality lets the microphones focus to the front for improved directional speech understanding especially in noisy environments – even for CIC wearers.



A world's first: Binaural OneMic directionality.

For enhanced natural directional hearing – even with CICs.

When technology improves nature: Sitting deep in the ear canal, CICs use the natural pinna effect for directional hearing to the front. Now equipped with the superior binax technology, even tiny hearing aids with only one microphone, such as CICs, can make use of e2e wireless 3.0 audio exchange technology. In a bilateral fitting, each hearing aid uses the information from its own microphone and from the other instrument's. This generates a virtual four-microphone network which allows for HDSR, and builds improved directionality to the front. This is what we call **binaural OneMic directionality**.

Adding to the natural pinna effect, binaural OneMic directionality provides the wearer with enhanced natural directional hearing.



All about binax.

If you would like to find out more about the groundbreaking binax technology, its features in detail, and the outstanding binax products, please take a look at our "binax" brochure.

The new Insio binax.

Directional hearing with the smallest hearing aids.

Insio binax introduces a world's first: Now hearing aid models with one microphone only can offer binaural directionality – even small CICs. And binax technology further enhances Insio binax: All models – from ITEs to tiniest IICs – are remotely controllable via the touchControl App.

For more details, please see page 26.



New!

The new Motion binax.

Sleek and easy-to-manage with true binaural capabilities.

The new Motion binax models cater to nearly any customer wish – and need. The rechargeable Motion SX binax and the DAI-ready Motion SA binax come in a new, compact housing design with improved directionality. And the rechargeable Motion PX binax provides a high level of amplification to cover a broad range of hearing losses.

For more details, please see page 24.



New!



For enhanced acoustic coupling, the new housings of Motion SX and SA binax have an interface with screw thread for the new ThinTube or redesigned earhook.

For more details, please see page 69.

New!

Unprecedented discretion:

Remote control option for smallest hearing aids.



touchControl App and Insio binax IIC

Equipped with the exceptional binax technology, now even our tiniest hearing aids, such as Ace binax™ and Insio binax IICs, are remotely controllable via the touchControl App.

New!


Performance level 3bx:

Now for all binax hearing instruments.

Offer the benefits of the outstanding binax features to almost all your customers: With Motion binax and Insio binax, we introduce the new binax performance level 3bx. And to help you further enhance your portfolio, we also present our first binax models Carat binax, Carat A binax, Pure binax™, and Ace binax in 3bx.

Hearing Solutions for Everyone.

We provide a broad range of hearing aids for every requirement and all degrees of hearing loss. You can offer your customers a wide selection from many different models and performance levels, from powerful RIC hearing aids, to versatile BTE hearing instruments or discreet ITE hearing aids. To facilitate orientation and help in the selection of the right hearing solution, all Siemens hearing instruments are placed in three categories according to their construction models. This way, you can quickly find out which model best suits a customer's wearing preferences.



RIC

With the receiver-in-canal, RIC models have a smaller housing that sits behind the ear. At the same time, they allow for sound transmission directly into the ear. This makes them ideal for wearers who would like a powerful and subtle solution with a high level of technological convenience.

BTE

Behind-the-ear models offer the optimum solution for customers who are at ease with their hearing aids and prefer a high level of wearing and handling comfort. We offer a large variety of BTEs for all levels of hearing loss.

ITE

For users who prefer a discreet solution, we offer a selection of tiny in-the-ear hearing aids that fit perfectly into the ear. These models are almost invisible and do not give the appearance that the customer is wearing a hearing instrument, but they still feature refined technology, power and convenience.

| Category | Hearing loss level | Available models | | | | | | | | | |
|----------|--------------------|------------------|---------|--|--|------------|------------|--|-------------|-------|--|
| | | RIC | profund | | | Pure binax | Pure micon | | Carat binax | Orion | |
| BTE | moderate to severe | | | | | | | | | | |
| | moderate | | | | | | | | | | |
| | mild | | | | | | | | | | |
| ITE | profund | | | | | | | | | | |
| | moderate to severe | | | | | | | | | | |
| | moderate | | | | | | | | | | |
| | mild | | | | | | | | | | |

Audio streaming and optional control app

Remotely controllable via app

Remotely controllable

Companion microphone option

Rechargeable























binax Feature Overview.

| Features | Information | Performance levels | | |
|---|---|--------------------|-------|-------|
| | | 7bx | 5bx | 3bx |
| General | Information | | | |
| Signal processing | Channels | 48 | 32 | 24 |
| Gain/MPO | Handles in Connexx | 20 | 16 | 12 |
| Hearing programs | Maximum available | 6 | 6 | 6 |
| touchControl App | iOS and Android devices only | ● | ● | ● |
| binaxFocus | | | | |
| Directionality | Resolution in channels | 48 | 32 | 24 |
| Binaural OneMic directionality ^B | Req. bilateral fitting & e2e wireless 3.0 | ●●●●● | — | — |
| Narrow directionality ^{A,C} | Req. bilateral fitting & e2e wireless 3.0 | ●●●●● | ●●●●○ | ●●●○○ |
| Directional microphones ^C | | ●●●●● | ●●●●○ | ●●●○○ |
| Spatial SpeechFocus ^{A,C} | Req. bilateral fitting & e2e wireless 3.0 | ●●●●● | — | — |
| SpeechFocus ^C | | ●●●●○ | ●●●●○ | — |
| TruEar™ | Req. BTE or RIC hearing aids | ●●●●● | ●●●●○ | ●●●○○ |
| Frequency compression | | ● | ● | ● |
| binaxSound | | | | |
| eWindScreen binaural ^{A,C} | Req. bilateral fitting & e2e wireless 3.0 | ● | — | — |
| eWindScreen | Adjustment steps | 3 | 3 | 1 |
| Extended bandwidth | | ● | — | — |
| Adaptive streaming volume ^A | Streaming only, req. easyTek & e2e wireless 3.0 | ● | — | — |
| SoundBrilliance™ ^A | Streaming only, req. easyTek & e2e wireless 3.0 | ● | ● | — |
| Feedback cancellation | | ● | ● | ● |
| Noise reduction | Resolution in channels | 48 | 32 | 24 |
| | Combined adjustment steps | 5 | 5 | 3 |
| Directional speech enhancement ^C | Individual adjustment steps | 3 | 1 | — |
| SoundSmoothing™ | Individual adjustment steps | 3 | 3 | 1 |
| Speech and noise management | Individual adjustment steps | 7 | 5 | 3 |
| binaxGuide | | | | |
| Sound equalizer | Number of classes | 6 | 3 | — |
| Data logging | | ● | ● | ● |
| Learning | Number of classes | 6 | 3 | 1 |
| Acclimatization manager | | ● | ● | ● |
| binax fit | | ● | ● | ● |
| Spatial Configurator ^C | Req. bilateral fitting & e2e wireless 3.0 | ● | ● | — |
| Span | Req. easyTek & easyTek App, or rocker switch | ● | ● | — |
| Direction | Req. easyTek & easyTek App | ● | ● | — |
| Tinnitus noiser | | | | |
| Static therapy signal (5 presets) | Handles in Connexx | 20 | 16 | 12 |
| Ocean wave therapy signal (4 presets) | | ● | ● | ● |

A) Not for Ace binax. B) Only for hearing aids with one microphone each. C) Only for hearing aids with two microphones each.









Advantages of a Higher Performance Level.

The following table gives you a short overview of the additional benefits that higher performance levels offer thanks to advanced features.

| Performance levels comparison | Listening situations where this feature may be helpful |
|--|---|
| 7bx benefits not available for 5bx | |
| 7bx' premium microphone technology in 48 channels with narrow directionality allows the wearer to understand speech even in difficult listening environments such as busy restaurants or cocktail parties. In some situations, they might even be able to hear better than a normal hearing person.* Even wearers of tiny Inzio binax single-microphone hearing aids can benefit from an improved directionality with the new binaural OneMic directionality . |   |
| With 7bx, the hearing instruments have special automatic settings for how wearers like to hear in the car. While suppressing car and traffic noise, spatial SpeechFocus allows them to hear what their fellow passengers are saying in the car, no matter if they are sitting in the front, back, left or right. |   |
| With 7bx' eWindScreen binaural , the wearer will have the highest listening comfort in windy situations. |  |
| No matter the background noise level, easyTek streaming remains loud and clear with 7bx' adaptive streaming volume . |   |
| 5bx benefits not available for 3bx | |
| 5bx' SpeechFocus makes it much easier for the wearer to follow a conversation in the car. |  |
| Thanks to 5bx's directional speech enhancement , the wearer may find it easier to follow conversations in noisy and crowded situations. |   |
| With 5bx, learning options enable the hearing aids to learn the wearer's listening preferences separately for three distinct acoustic situations: Speech, noise, and music. |   |
| 3bx benefits not available for Orion | |
| 3bx Bluetooth compatibility allows the wearer to stream phone conversations, TV, and music directly and in true stereo into their ears, without annoying interference or whistling. |   |
| With 3bx' 24 channel directional microphone , including narrow directionality , and noise reduction processing, background noises are effectively reduced. |  |
| 3bx' new ocean wave sounds expand the flexibility of the tinnitus therapy solution and create a more positive, soothing, and stress-relieving listening experience for people with tinnitus. |  |
| Orion benefits not available for Sirion | |
| Depending on the wearer's hearing loss, Orion's frequency compression allows the wearer to hear high-pitched sounds such as consonant speech sounds and phone rings, that were not possible before. |   |
| Orion's superior fitting scheme offers the wearer a more natural sound impression. |   |
| Sirion benefits not available for Basic | |
| Sirion's 6 frequency shaping and compression channels offer more fitting flexibility for the wearer's hearing loss. |  |
| With Sirion's feedback cancellation , the wearer may encounter less annoying whistling. |  |

* See legal note on final page in this document.

Listening situations improved by advanced features:

-  At home
-  Outdoor conversation
-  Restaurant
-  Phone calls/watching TV
-  Car
-  Social events
-  Meetings and conferences
-  Music

miCON, Essential and Basic Feature Overview.

| Features | Performance levels | | | | | | |
|--|--------------------|---------|---------|----------------|----------------|-----------|-------|
| | 7mi | 5mi | 3mi | Essential | | Basic | |
| | | | | Orion | Sirion | Intuis | Lotus |
| Better hearing | | | | | | | |
| Frequency channels/Compression controls | 48/20 ^B | 32/16 | 24/12 | 16/8 | 12/6 | 4/4 | 2/2 |
| Number of programs | 6 | 6 | 6 | 4 | 4 | 1-4 | 1-3 |
| miFocus™ A | ● | ● | ● | - | - | - | - |
| High res. adaptive TwinMic System™ | 48 | 32 | 24 | 16 | - | - | - |
| Directional speech enhancement | ● | - | - | - | - | - | - |
| Frequency compression | ● | ● | ● | ● | - | - | - |
| High res. SpeechFocus | ● | ● | - | - | - | - | - |
| micon TruEar | ● | ● | - | - | - | - | - |
| Directional microphones ^A | - | - | - | Automatic | ● ^C | Automatic | - |
| Sound comfort | | | | | | | |
| miSound™ | ● | ● | ● | - | - | - | - |
| Extended bandwidth (up to 12 kHz) | ● | - | - | - | - | - | - |
| Feedback cancellation | ● | ● | ● | ● | ● | - | - |
| micon speech and noise management | 7 steps | 5 steps | 3 steps | - | - | - | - |
| micon SoundSmoothing | 3 steps | 3 steps | on/off | - | - | - | - |
| micon eWindScreen | 3 steps | 3 steps | on/off | - | - | - | - |
| micon SoundBrilliance | 3 steps | 3 steps | - | - | - | - | - |
| FeedbackBlocker™ | - | - | - | - | - | ● | ● |
| Speech and noise management | - | - | - | 3 steps | on/off | on/off | - |
| eWindScreen | - | - | - | ● | - | - | - |
| e2e wireless 2.0 | ● | ● | ● | ● ^D | - | - | - |
| Fitting/Individuality | | | | | | | |
| miGuide™ | ● | ● | ● | - | - | - | - |
| Acclimatization manager | ● | ● | - | - | - | - | - |
| micon learning | ● | ● | ● | - | - | - | - |
| Sound equalizer/no. of situations | 6 | 3 | 1 | - | - | - | - |
| micon fit | ● | ● | ● | - | - | - | - |
| Data logging | ● | ● | ● | ● | ● | - | - |
| Audio streamer/Instrument control | | | | | | | |
| easyPocket | ● | ● | ● | ● | - | - | - |
| miniTek and miniTek Remote App | ● | ● | ● | - | - | - | - |
| ConnexLink | ● | ● | ● | ● | ● | - | - |

A) Not for CIC. B) 24/16 for Nitro micon BTE. C) Only BTE. D) Left and right synchronization.

Advantages of a Higher Performance Level.

The following table gives you a short overview of the additional benefits that higher performance levels offer thanks to advanced features.

| Performance levels comparison | Listening situations where this feature may be helpful |
|---|--|
| 7mi benefits not available for 5mi | |
| Depending on their hearing loss, the wearer may appreciate a naturally brilliant sound quality with 7mi's extended bandwidth . | |
| Thanks to 7mi's directional speech enhancement , the wearer may find it easier to follow conversations in noisy and crowded situations. | |
| With miGuide that detects and differentiates 6 distinct sound environments, 7mi hearing instruments offer optimal settings for each of these environments. Furthermore, 7mi instruments can learn wearer listening preferences for these different environments. | |
| 5mi benefits not available for 3mi | |
| Thanks to 5mi's SpeechFocus , it is much easier for the wearer to follow a conversation in the car. | |
| With miGuide that detects and differentiates 3 distinct sound environments, 5mi hearing instruments offer optimal settings for each of these environments. Furthermore, 5mi can learn wearer listening preferences for these different environments. | |
| 3mi benefits not available for Orion | |
| 3mi's Bluetooth compatibility allows the wearer to stream phone conversations, TV, and music directly and in true stereo into their ears, without annoying interference or whistling. | |
| Thanks to 3mi's eWindScreen , the wearer may have an easier time following conversations in windy outdoor situations. | |
| Orion benefits not available for Sirion | |
| Depending on the wearer's hearing loss, Orion's frequency compression allows the wearer to hear high-pitched sounds such as consonant speech sounds and phone rings, that were not possible before. | |
| Orion's superior fitting scheme offers the wearer a more natural sound impression. | |
| Sirion benefits not available for Basic | |
| Sirion's 6 frequency shaping and compression channels offer more fitting flexibility for the wearer's hearing loss. | |
| With Sirion's feedback cancellation , the wearer may encounter less annoying whistling. | |

Listening situations improved by advanced features:

- Conversations
- Meetings and conferences
- Music
- At home
- Outdoor conversation
- Restaurant
- Phone calls/watching TV
- Car
- Social events

Ace binax.

A highlight of discretion and sound.

Key Features

- Small size
- Remotely controllable
- Battery size 10A
- Optimized ergonomic design for premium performance

As our tiniest RIC, Ace has always set a benchmark in minimizing size and maximizing discretion. Still so small that it is almost invisible when worn, the next generation of BestSound Technology, further enhances Ace binax' level of discretion. Remotely controllable via the touchControl App, your customers can operate their hearing aid programs, volume, bass and treble with a few quick taps on their smartphones. This makes Ace binax an optimum solution for first-time wearers or people that demand style and convenience first. With Ace binax, your customers do not have to choose between discretion and functionality – now they can have both.

Ace binax. Now more discreet than ever.

Remotely controllable
Via the touchControl App, hearing programs, volume, bass and treble can be discreetly changed.

Exchangeable housing
Housing can be easily exchanged.

Ultra small design
The hearing instrument is nearly imperceptible when worn.


Directional microphones
Optimized placement and orientation ensure the best possible hearing performance.

IP67-rated
Resistant to moisture, sweat, dust and dirt.

Push button
Also available as flat cover.

New! Tinnitus noiser
Tinnitus therapy signals only or mixed mode. With five static noise types and four ocean wave signals.

miniReceiver 2.0
For excellent sound performance. Available in four different lengths and power levels S, M, P and HP.



Accessories, see pages 56–63. Accessory compatibility overview, see page 57.



Colors:




Ace binax

| Type | RIC | | | |
|------------------------------|-------------------------|--------------------|--------------------|---------------------|
| Performance levels | 7bx/5bx/3bx New! | | | |
| Description | miniReceiver 2.0 S | miniReceiver 2.0 M | miniReceiver 2.0 P | miniReceiver 2.0 HP |
| Picture | | | | |
| Fitting range | | | | |
| Product details | | | | |
| Matrix (2 ccm) | RIC 45/108 dB | RIC 55/113 dB | RIC 60/118 dB | RIC 65/123 dB |
| Battery size | 10A | 10A | 10A | 10A |
| Power consumption | 0.9 mA | 1.0 mA | 0.9 mA | 0.9 mA |
| Features | | | | |
| BestSound Technology | binax | binax | binax | binax |
| e2e wireless 3.0 | - | - | - | - |
| touchControl App | ● | ● | ● | ● |
| Rechargeable | - | - | - | - |
| Volume control | - | - | - | - |
| Push button | ● ⁽¹⁾ | ● ⁽¹⁾ | ● ⁽¹⁾ | ● ⁽¹⁾ |
| T-coil | - | - | - | - |
| Battery door on/off function | ● | ● | ● | ● |
| Audio input | - | - | - | - |
| Tinnitus noiser | ● | ● | ● | ● |
| Modifications | | | | |
| Color conversion kit | ○ | ○ | ○ | ○ |
| Accessories | | | | |
| Programming adapter | size 10 | size 10 | size 10 | size 10 |
| miniReceiver fitting set | 2.0 | 2.0 | 2.0 | 2.0 |

(1) Exchangeable controls – programmable via Connexx

● Standard
○ Optional




miniReceiver 2.0 measuring gauge
With the miniReceiver 2.0 measuring gauge, you can easily determine whether your customer can be fitted with a power miniReceiver 2.0 HP. See details on page 67.

Pure binax.

Sophisticated binaural technology in an elegant design.

- Key Features**
- Rechargeable
 - Fully featured
 - Battery size 312
 - Wireless connectivity


So small, yet so fully-featured: Equipped with binax, the next generation of BestSound Technology, Pure binax is ready to delight even the most discerning wearers. It offers all the outstanding advantages of a binaural hearing instrument whenever needed because it belongs to the industry's most energy efficient hearing aids when it comes to binaural audio transmission. Plus, two independent clinical studies have proven that Pure binax allows wearers to hear in difficult listening situations such as restaurants or parties even better than people with normal hearing.*
Pure binax. Discretion meets a new dimension of hearing.



Pure binax

| Type | RIC | | | |
|------------------------------|-------------------------|--------------------|--------------------|---------------------|
| Performance levels | 7bx/5bx/3bx New! | | | |
| Description | miniReceiver 2.0 S | miniReceiver 2.0 M | miniReceiver 2.0 P | miniReceiver 2.0 HP |
| Picture | | | | |
| Fitting range | | | | |
| Product details | | | | |
| Matrix (2 ccm) | RIC 45/108 dB | RIC 60/119 dB | RIC 70/124 dB | RIC 75/130 dB |
| Battery size | 312 | 312 | 312 | 312 |
| Power consumption | 0.9 mA | 1.0 mA | 1.0 mA | 1.1 mA |
| Features | | | | |
| BestSound Technology | binax | binax | binax | binax |
| e2e wireless 3.0 | ● | ● | ● | ● |
| touchControl App | ● | ● | ● | ● |
| Rechargeable | ● | ○ ² | ○ ² | ○ ² |
| Volume control | ○ ¹ | ○ ¹ | ○ ¹ | ○ ¹ |
| Push button | ○ ¹ | ○ ¹ | ○ ¹ | ○ ¹ |
| T-coil | ● | ● | ● | ● |
| Battery door on/off function | ● | ● | ● | ● |
| Audio input | ○ ³ | ○ ³ | ○ ³ | ○ ³ |
| Tinnitus noiser | ● | ● | ● | ● |
| Modifications | | | | |
| Color conversion kit | ● | ● | ● | ● |
| Accessories | | | | |
| Programming adapter | size 312 | size 312 | size 312 | size 312 |
| miniReceiver fitting set | 2.0 | 2.0 | 2.0 | 2.0 |

(1) Exchangeable controls – programmable via Connexx (2) Reduced battery life (3) With easyTek
 ● Standard ○ Optional



* See legal note on final page in this document.

Carat binax.

True binaural hearing meets flexibility in a compact design.

Key Features

- Rechargeable
- Option: DAI-ready
- Battery size 13
- Wireless connectivity

Carat binax™ delivers a natural binaural listening experience by combining balanced design with binax technology. Two independent clinical studies have proven that Carat binax allows wearers to hear even better than people with normal hearing in certain listening situations.* Both the rechargeable Carat binax and Carat A binax with optional DAI offer binax' binaural benefits and wireless connectivity, flexibly meeting customer needs.

Carat binax. Flexibility meets a new dimension of hearing.

e2e wireless 3.0
This third-generation wireless data exchange system is what enables the sophisticated binaural features in binax.

Directional microphones
Optimized placement and orientation ensure the best possible hearing performance.

Exchangeable housing
Housing can be easily exchanged.

IP67-rated
Resistant to moisture, sweat, dust and dirt.

Rocker switch
For program and volume adjustments.

Integrated T-coil
For inductive transmission.

Tinnitus noiser **New!**
Tinnitus noiser function only or mixed mode. With five static noise types and four ocean wave signals.

Wireless connectivity and remotely controllable
Easily operable via the easyTek audio streamer and remote control plus the touchControl App.

Charging contact
No hassle with batteries. The eCharger charges and dries Carat binax overnight.

miniReceiver 2.0
For excellent sound performance. Available in four different lengths and power levels S, M, P and HP.

Flexible battery compartment
For 13 standard battery, Carat binax models are rechargeable. Carat A binax models come with a standard battery compartment that can be exchanged for one with a connection for direct audio input.

Carat binax

Carat A binax


Actual size in cm

Carat binax

Carat A binax with audio shoe

* See legal note on final page in this document.


Carat binax Carat A binax



| Type | RIC | | | |
|------------------------------|-------------------------|--------------------|--------------------|---------------------|
| Performance levels | 7bx/5bx/3bx New! | | | |
| Description | miniReceiver 2.0 S | miniReceiver 2.0 M | miniReceiver 2.0 P | miniReceiver 2.0 HP |
| Picture | | | | |
| Fitting range | | | | |
| Product details | | | | |
| Matrix (2 ccm) | RIC 45/108 dB | RIC 60/119 dB | RIC 70/124 dB | RIC 75/130 dB |
| Battery size | 13 | 13 | 13 | 13 |
| Power consumption | 0.9 mA | 1.1 mA | 1.0 mA | 1.1 mA |
| Features | | | | |
| BestSound Technology | binax | binax | binax | binax |
| e2e wireless 3.0 | ● | ● | ● | ● |
| touchControl App | ● | ● | ● | ● |
| Rechargeable ⁽⁴⁾ | ● | ● | ● | ○ ⁽⁵⁾ |
| Volume control | ● | ● | ● | ● |
| Push button | - | - | - | - |
| T-coil | ● | ● | ● | ● |
| Battery door on/off function | ● | ● | ● | ● |
| Audio input ⁽⁵⁾ | ○ ⁽³⁾ | ○ ⁽³⁾ | ○ ⁽³⁾ | ○ ⁽³⁾ |
| Tinnitus noiser | ● | ● | ● | ● |
| Modifications | | | | |
| Color conversion kit | ● | ● | ● | ● |
| Accessories | | | | |
| Programming adapter | size 13 | size 13 | size 13 | size 13 |
| miniReceiver fitting set | 2.0 | 2.0 | 2.0 | 2.0 |

(2) Reduced battery life
(3) With easyTek
(4) Carat binax models only
(5) Carat A binax models only

● Standard
○ Optional



miniReceiver 2.0 measuring gauge
With the miniReceiver 2.0 measuring gauge, you can easily determine whether your customer can be fitted with a power miniReceiver 2.0 HP. See details on page 67.

Motion binax.

New!

Sleek and easy to manage with true binaural capabilities to fit any ear.

Key Features

- Sleek and compact housing
- Easy handling
- Fully featured

The new Motion binax models offer true binaural capabilities. Fully featured, they provide almost all options to cover any individual need. The rechargeable Motion SX binax and the DAI-ready Motion SA binax come in a new, sleek design with improved directionality. This makes them not only very discreet and comfortable to wear but also provides the best possible listening performance. The new Motion PX binax is a powerful partner in a proven design: Due to its high level of amplification, it can cover a broad range of hearing losses and is a reliable solution for almost any customer. Its user-friendly controls and intelligent interfaces offer wearers hassle-free handling. Thanks to their outstanding wireless connectivity, all Motion binax models offer your customers a maximum of convenience. **Motion binax. Comfort meets a new dimension of hearing.**

Directional microphones
Optimized placement and orientation provide an improved directionality to ensure the best possible binaural listening performance.

Wireless connectivity
Easy audio streaming and remote control functionality via the easyTek and the touchControl App.

New! Tinnitus noiser
Tinnitus therapy signals only or mixed mode. With five static noise types and four ocean wave signals.

e2e wireless 3.0
This third-generation wireless data exchange system is what enables the sophisticated binaural features in binax.

Charging contact
The eCharger charges and dehumidifies Motion SX and PX binax models overnight.

ThinTube New!
The ultra-slim, ultra-transparent ThinTube comes with a new concha lock position. It is unobtrusive, easy to insert, and combines stability with absolute flexibility.

Option: Direct audio input
Motion SA binax models come with a standard battery compartment that can be exchanged for one with a direct audio input connection.

Motion SX binax

Motion SA binax

Motion SX/SA binax

Motion PX binax

Actual size in cm

Motion binax

Colors:



Accessories, see pages 56–63. Accessory compatibility overview, see page 57.



Additional colors SX and SA



Additional colors SA



Additional colors PX



Motion binax

| Type | BTE | | |
|------------------------------|---|---|---|
| Performance levels | 7bx/5bx/3bx | | |
| Description | Motion SX binax | Motion SA binax | Motion PX binax |
| Picture | | | |
| Fitting range | | | |
| Product details | ThinTube 53/125 dB Earhook 60/130 dB | ThinTube 53/125 dB Earhook 60/130 dB | S-LifeTube 65/127 dB Earhook 75/135 dB |
| Matrix (2 ccm) | | | |
| Battery size | 13 | 13 | 13 |
| Power consumption | 1.0 mA | 1.0 mA | 1.1 – 1.3 mA |
| Features | | | |
| BestSound Technology | binax | binax | binax |
| e2e wireless 3.0 | ● | ● | ● |
| touchControl App | ● | ● | ● |
| Rechargeable | ● | - | ● |
| Volume control | ● ⁷ | ● ⁷ | ● |
| Push button | - | - | ● |
| T-coil | ● | ● | ● |
| Battery door on/off function | ● | ● | ● |
| Battery door lock | - | - | - |
| Audio input | ○ ³ | ○ ³ | ○ ³ |
| Tinnitus noiser | ● | ● | ● |
| Modifications | | | |
| Color conversion kit | ● | ● | ● |
| Eyeglass adapter | - | - | ○ |
| Accessories | | | |
| Programming adapter | size 13 | size 13 | size 13 |
| Fitting set | ThinTube | ThinTube | S-LifeTube |

(3) With easyTek
(7) Programmable via Connexx
(8) With battery compartment modification

● Standard
○ Optional

Insió binax.

New!

Leading in binaural directionality.

Key Features

- Binaural directionality even in CICs
- Small construction due to proprietary components
- Remote control for all models

The new Insió binax presents an astonishing world's first: Thanks to the superior binax technology, now even hearing aids with one microphone can offer true binaural directionality. Sitting in the ear canal, standard CICs benefit from the directional effect of the outer ear. Now thanks to binaural OneMic directionality, Insió binax CICs can even improve this directional effect. The result is a directionality that cannot be achieved without binaural technology. Additionally, all Insió models, even tiniest IICs, are remotely controllable via the touchControl App. This way, the new Insió binax can offer a nearly invisible solution with outstanding audiological performance.

Insió binax. Minimized design, maximized performance.

New! Wireless connectivity
This third-generation wireless data exchange system is what enables the sophisticated binaural features in binax.

Vent
For optimum ear ventilation. With Optivent option.


Battery door

Remotely controllable New!
The hearing programs, volume, bass and treble of all Insió binax models can be discreetly changed via the touchControl App.

Tinnitus noiser New!
Tinnitus therapy signals only or mixed mode. With five static noise types and four ocean wave signals.

Directional microphones
Optimized placement and orientation ensure the best possible binaural hearing performance.

New! Binaural OneMic directionality
In a bilateral fitting, even hearing aids with only one microphone can offer enhanced natural directional hearing.



Actual size in cm



Accessories, see pages 56–61. Accessory compatibility overview, see page 57.

Colors IIC:
Dark brown faceplate with white shell.

Standard Skin Colors CIC, ITC, ITE:
See details on page 91.

Insió binax

| Type | ITE | | | |
|----------------------|-------------------------|-------------------------|-------------------------|-------------------------|
| Performance levels | 7bx/5bx/3bx | | | |
| Description | IIC | CIC | ITC | ITE |
| Picture | | | | |
| Fitting range | | | | |
| Product details | | | | |
| Matrix (2 ccm) | 50/113 ¹¹ dB | 50/113 ¹¹ dB | 55/118 ¹¹ dB | 60/123 ¹¹ dB |
| Battery size | 10 | 10 | 312/10 | 13/312 |
| Power consumption | 0.9 mA | 1.0 mA | 1.1 mA | 1.1 mA |
| Features | | | | |
| BestSound Technology | binax | binax | binax | binax |
| e2e wireless 3.0 | - | ○ | ○ | ○ |
| touchControl App | ● | ● | ● | ● |
| Volume control | - | ○ | ○ | ○ |
| Push button | - | ○ | ○ | ○ |
| T-coil | - | - | ○ | ○ |
| Audio input | - | ○ ³ | ○ ³ | ○ ³ |
| Tinnitus noiser | - | ● | ● | ● |
| Optivent | ● | ● | ● | ● |
| Accessories | | | | |
| Programming adapter | Flex connector | Flex connector | Flex connector | Flex connector |

(3) With easyTek
(11) Max. values

● Standard
○ Optional

Pure micon.

Discreet and versatile.

Key Features


- Rechargeable
- Fully featured
- Discreet design
- Wireless connectivity

With a very discreet, understated look, this rechargeable and wireless solution comes in a variety of models. Pure micon™ – our fully featured RIC, has an exceptionally functional design, offering the best audiological performance in a very small size. Pure micon models cover practically all levels of hearing loss, and deliver outstanding versatility in a RIC product.

Pure micon. Perfectly equipped for effortless interaction.



For more information, please consult our "Pure micon" brochure.




Colors:



Pure micon

| Type | RIC | | | |
|------------------------------|--------------------|--------------------|--------------------|---------------------|
| Performance levels | 7mi/5mi/3mi | | | |
| Description | miniReceiver 2.0 S | miniReceiver 2.0 M | miniReceiver 2.0 P | miniReceiver 2.0 HP |
| Picture | | | | |
| Fitting range | | | | |
| Product details | | | | |
| Matrix (2 ccm) | RIC 45/108 dB | RIC 60/119 dB | RIC 70/124 dB | RIC 75/130 dB |
| Battery size | 312 | 312 | 312 | 312 |
| Power consumption | 0.8 mA | 1.0 mA | 1.0 mA | 1.1 mA |
| Features | | | | |
| BestSound Technology | micon | micon | micon | micon |
| e2e wireless 2.0 | ● | ● | ● | ● |
| Rechargeable | ● | ○ ² | ○ ² | ○ ² |
| Volume control | ○ ¹ | ○ ¹ | ○ ¹ | ○ ¹ |
| Push button | ○ ¹ | ○ ¹ | ○ ¹ | ○ ¹ |
| T-coil | ● | ● | ● | ● |
| Battery door on/off function | ● | ● | ● | ● |
| Audio input | ○ ⁶ | ○ ⁶ | ○ ⁶ | ○ ⁶ |
| Tinnitus noiser | ● | ● | ● | ● |
| Modifications | | | | |
| Color conversion kit | ● | ● | ● | ● |
| Accessories | | | | |
| Programming adapter | size 312 | size 312 | size 312 | size 312 |
| miniReceiver fitting set | 2.0 | 2.0 | 2.0 | 2.0 |

- (1) Exchangeable controls – programmable via Connexx
 (2) Reduced battery life With miniTek
 ● Standard
 ○ Optional



miniReceiver 2.0 measuring gauge
 With the miniReceiver 2.0 measuring gauge, you can easily determine whether your customer can be fitted with a power miniReceiver 2.0 HP. See details on page 67.


Siemens Life micon.

Sense of hearing meets sense of style.

Key Features

- Fully automatic functions
- Tinnitus option
- Wireless connectivity

Siemens Life micon™ is one of the smallest standard BTE in the Siemens product range and the ideal alternative if a RIC hearing aid is not an option. It is appropriate for a wide variety of ear canals, whether bent or simply extremely small, and also if a customer has cerumen issues. Fitting is simple thanks to the flexible LifeTubes – Siemens Life micon fits virtually every user right away. It is simple, convenient and reliable – also thanks to its proven housing. **Siemens Life micon. The small work of art that's big on style.**



For more information, please consult our "Siemens Life micon" brochure.



Accessories, see pages 56–63.
Accessory compatibility overview, see page 57.



Colors:



Siemens Life micon

| Type | BTE | |
|------------------------------|--------------------|-------------------|
| Performance levels | 7mi/5mi/3mi | |
| Description | Siemens Life micon | |
| Picture | | |
| Fitting range | | |
| Product details | | |
| Matrix (2 ccm) | LifeTube 47/121 dB | Earhook 57/123 dB |
| Battery size | 312 | 312 |
| Power consumption | 0.8 mA | 0.8 mA |
| Features | | |
| BestSound Technology | micon | micon |
| e2e wireless 2.0 | ● | ● |
| Rechargeable | - | - |
| Volume control | ○ ⁽⁷⁾ | ○ ⁽⁷⁾ |
| Push button | ● | ● |
| T-coil | ○ ⁽⁶⁾ | ○ ⁽⁶⁾ |
| Battery door on/off function | ● | ● |
| LifeTube/Standard earhook | ● | ○ |
| Audio input | ○ ⁽⁶⁾ | ○ ⁽⁶⁾ |
| Tinnitus noiser | ● | ● |
| Modifications | | |
| Color conversion kit | ● | ● |
| Accessories | | |
| Programming adapter | size 312 | size 312 |
| LifeTube fitting set | ● | ● |

(6) With miniTek
(7) Programmable via Connex

● Standard
○ Optional


Motion micon.

A partner for life.

Key Features

- Onboard controls
- Wireless connectivity
- Fully featured to suit all ages

Motion micon™ can fit any ear and is well-suited for active lifestyles: Fully featured, hassle-free, reliable and with a high degree of wearing and handling comfort. This family comes in four versions, which cover almost all degrees of hearing loss and all age groups. The M and P models include many safety features perfect for infants and toddlers. The SX and PX versions offer the convenience of rechargeability. And all Motion micon models are IP67-rated. **Motion micon. Superb sound meets outstanding comfort.**



For more information, please consult our "Motion micon" brochure.

Slim-neck design
Discreet fit behind the ear, high degree of wearing comfort.

Directional microphones
Optimized placement and orientation ensure the best possible hearing performance.

Dual onboard controls
Push button and rocker switch on one hearing aid.

Integrated T-coil
For inductive transmission.

Tinnitus noiser
Two programmable modes: Tinnitus therapy function only or mixed mode.

IP67-rated
Resistant to moisture, sweat, dust and dirt.

Superior nanocoating technology
Protects the hearing aid from corrosion and water damage.

Charging contact
No hassle with batteries. Available for SX and PX models.


Option: Direct audio input
Motion M and P micon models are equipped with a child-lock battery door so children don't unintentionally open it. Modified, these models easily connect to most common FM solutions via an optional audio shoe.

Motion PX micon




Motion M/P micon

Actual size in cm



Motion micon

Colors:



Accessories, see pages 56–63.
Accessory compatibility overview, see page 57.



Motion micon

| Type | BTE | | | |
|------------------------------|---|---|---|---|
| Performance levels | 7mi/5mi/3mi | | | |
| Description | Motion SX micon | Motion M micon | Motion P micon | Motion PX micon |
| Picture | | | | |
| Fitting range | | | | |
| Product details | | | | |
| Matrix (2 ccm) | LifeTube 53/125 dB Earhook 60/130 dB | S-LifeTube 60/124 dB Earhook 65/130 dB | S-LifeTube 65/127 dB Earhook 75/135 dB | S-LifeTube 65/127 dB Earhook 75/135 dB |
| Battery size | 13 | 13 | 13 | 13 |
| Power consumption | 1.0 mA | 1.1 mA | 1.1–1.3 mA | 1.1–1.3 mA |
| Features | | | | |
| BestSound Technology | micon | micon | micon | micon |
| e2e wireless 2.0 | ● | ● | ● | ● |
| Rechargeable | ● | - | - | ● |
| Volume control | ○ ¹ | ● | ● | ● |
| Push button | ○ ¹ | ● | ● | ● |
| T-coil | ● | ● | ● | ● |
| Battery door on/off function | ● | ● | ● | ● |
| LifeTube/Standard earhook | ● | ● | ● | ● |
| Battery door lock | - | ○ | ○ | - |
| Audio input | ○ ⁶ | ○ ⁶ | ○ ⁶ | ○ ⁶ |
| Tinnitus noiser | ● | ● | ● | ● |
| Modifications | | | | |
| Color conversion kit | ● | ● | ● | ● |
| Eyeglass adapter | ○ | ○ | ○ | ○ |
| Accessories | | | | |
| Programming adapter | size 13 | size 13 | size 13 | size 13 |
| Fitting set | LifeTube | S-LifeTube | S-LifeTube | S-LifeTube |

(1) Exchangeable controls – programmable via Connexx
(6) With miniTek
(8) With battery compartment modification

● Standard
○ Optional

Aquaris micon.


Waterproof, robust and more powerful than ever.

Showering, gardening, biking, walking in the rain or swimming – Aquaris micon™ can handle it all. The housing of Aquaris micon is completely protected, making it suitable for continuous immersion in water – the first hearing aid of its kind to be successfully tested to the IP68-standard. With 65 dB of gain and micon generation technology, Aquaris micon offers more power than ever. Sweat, dust, dirt and sand mean nothing to the tough and durable Aquaris micon – perfect for anyone who wants to enjoy life without limits.

Aquaris micon. Lets sound in, not the elements.

Key Features

- Water- and dustproof
- Robust construction
- Secure fit



For more information, please consult our "Aquaris micon" brochure.



GORE™ Acoustic Vents for Microphone Protection*

The high-performance acoustic membrane preserves the microphone's sound quality while protecting it from water and dust particles.

Waterproof microphone membrane

The microphone membranes are both waterproof and "sound-permeable", enabling the highest sound quality possible.

Microphone cover

Not removable, it protects the microphone membrane and reduces unpleasant noises, such as wind, raindrops, or the rustling of hair.

Improved battery door

The durable silicone-sealed battery door is improved to prevent direct membrane contact when opening the battery compartment.

Receiver protection

A waterproof gasket at the housing opening stops water drops from entering.

Tinnitus noiser

Two programmable modes: Tinnitus therapy function only or mixed mode.

Bottle housing

The exchangeable housing shells are produced as one complete piece and connected to the battery compartment using a silicone sealant.

IP68-rated

Completely waterproof and dustproof.

Soft-touch surface


For a better grip behind the ear with anti-repelling effect.

Improved waterproof membrane for air exchange

The larger membrane ensures even better airflow to the battery while its smooth surface is designed to reliably repel water.



Aquaris micon



* GORE and designs are trademarks of W. L. Gore & Associates.

Accessories, see pages 56 – 63. Accessory compatibility overview, see page 57.



Colors:



Aquaris micon

| Type | BTE | |
|------------------------------|--------------------|-------------------|
| Performance levels | 7mi/5mi/3mi | |
| Description | Aquaris micon | |
| Picture | | |
| Fitting range | | |
| Product details | | |
| Matrix (2 ccm) | LifeTube 49/122 dB | Earhook 65/129 dB |
| Battery size | 13 | 13 |
| Power consumption | 1.1 mA | 1.0 mA |
| Features | | |
| BestSound Technology | micon | micon |
| e2e wireless 2.0 | ● | ● |
| Rechargeable | - | - |
| Volume control | ○ ⁹ | ○ ⁹ |
| Push button | ○ ⁹ | ○ ⁹ |
| T-coil | ○ ⁶ | ○ ⁶ |
| Battery door on/off function | ● | ● |
| LifeTube/Standard earhook | ● | ● |
| Audio input | ○ ⁶ | ○ ⁶ |
| Tinnitus noiser | ● | ● |
| Modifications | | |
| Color conversion kit | ● | ● |
| Accessories | | |
| Programming adapter | CS 44 | CS 44 |
| Fitting set | ● | ● |

(6) With miniTek ● Standard
(9) With remote control ○ Optional



1) aquapac

aquapac is the ideal place to store remote controls, mobile phones and mp3 players when at the pool or in the garden. With up to 1 m water pressure resistance, it ensures devices are well protected. Better still, mp3 players can be placed in the aquapac and used while swimming.

Item No. 10668262

2) Sport Clip

The Sport Clip ensures a secure hold in any situation.

Item No. 10666837

Nitro micon.


Maximum power for maximum gain.

Key Features

- Maximum amplification, even in the directional microphone mode
- Wireless connectivity

Nitro micon™ is the ideal solution for people with severe to profound hearing loss: It offers more speech understanding, sound quality and wearing comfort to those who need it most. In addition, nearly all models offer the unique combination of maximum amplification with wireless connectivity. By combining super power with sophisticated technology, Nitro micon provides the options, functions and amplification that wearers who are dependent on their hearing aids need to stay connected to their world.

Nitro micon. With super power and super connectivity.



For more information, please consult our "Nitro micon" brochure.

Directional microphones
Optimized placement and orientation ensure the best possible directivity.

Ergonomic slim-neck design
Discreet fit behind the ear, high degree of wearing comfort.

Dual onboard controls
Integrated push button and rocker switch for manual control including a programmable key lock to prevent unintentional program or volume changes.

AutoPhone
Automatically switches to telephone mode when the telephone is held to the ear.

IP67-rated
Resistant to moisture, sweat, dust and dirt.

LED status indicator
Shows instrument's program change, key lock activation and deactivation.


Option: Direct audio input
Nitro micon is equipped with a child-lock battery door so children don't unintentionally open it. Modified, these models easily connect to most common FM solutions via an optional audio shoe.

Tinnitus noiser
Two programmable modes: Tinnitus therapy function only or mixed mode.

Directional microphones*
Focus on sound from in front of the wearer while reducing noises from other directions.

Vent
For optimum ear ventilation. With Optivent option.

Nitro micon BTE Nitro micon ITE



Frequency channels
The input is received and processed in 48 channels, but due to receiver restrictions only 24 channels are available at the output stage. Frequency compression can be used to make the whole input range audible.

* Available for most models.

Accessories, see pages 56–63. Accessory compatibility overview, see page 57.



Colors:

BTE



ITE




Nitro micon

| Type | BTE | ITE |
|------------------------------|-------------------|---|
| Performance levels | 7mi/3mi | 7mi/3mi |
| Description | Nitro micon | Nitro micon CIC, ITC, ITE |
| Picture | | |
| Fitting range | | |
| Product details | | |
| Matrix (2 ccm) | Earhook 82/141 dB | for all models from 55/118 dB ¹¹ up to 71/130 ¹¹ dB |
| Battery size | 675 | CIC 10, ITC 312/10, ITE 13/312 |
| Power consumption | 2.1 mA | 1.0 – 1.2 mA |
| Features | | |
| BestSound Technology | micon | micon |
| LED status indicator | ● ⁷ | - |
| e2e wireless 2.0 | ● | ○ |
| Rechargeable | - | - |
| Volume control | ● | ○ |
| Push button | ● ⁷ | ○ |
| T-coil | ● | ○ ¹⁰ |
| Battery door on/off function | ● | ● |
| Battery door lock | ○ ⁸ | - |
| Audio input | ○ ⁸ | ○ ⁶ |
| Tinnitus noiser | ● | ● |
| Optivent | - | ○ |
| Modifications | | |
| Color conversion kit | only top housing | - |
| Accessories | | |
| Programming adapter | CS 44 | Flex connector |

(6) With miniTek
 (7) Programmable via Connexx
 (8) With battery compartment modification
 (10) Not available for CIC
 (11) Max. values

● Standard
 ○ Optional



Nitro micon measuring gauge
The measuring gauge supports the fitting process by letting you easily find out about the right size for the wearer's custom instruments. See details on page 67.

Insió micon.

Engineered for absolute comfort and discretion.

Key Features

- Optivent option
- Wireless connectivity




For more information, please consult our "Insió micon" brochure.

Our ultra-tiny, ultra-comfortable Insió micon™ is a discreet solution with micon technology – optimized to fit the specific characteristics of ITE hearing aids. It offers an individualized hearing experience and Optivent™ ventilation to support a perfect balance between speech understanding and sound quality. Insió micon's technology and design provide a high level of discretion and wearing comfort, outstanding listening satisfaction and increased understanding from the start.

Insió micon. Live your individuality.



Actual size in cm



Insió micon



Accessories, see pages 56 – 63. Accessory compatibility overview, see page 57.

Colors IIC:
Dark brown faceplate with white shell.

Standard Skin Colors CIC, ITC, ITE:
See details on page 91.

| Type | ITE | | | |
|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|
| Performance levels | 7 mi/5 mi/3 mi | | | |
| Description | IIC | CIC | ITC | ITE |
| Picture | | | | |
| Fitting range | | | | |
| Product details | | | | |
| Matrix (2 ccm) | 50/113 ⁽¹¹⁾ dB | 50/113 ⁽¹¹⁾ dB | 55/118 ⁽¹¹⁾ dB | 60/123 ⁽¹¹⁾ dB |
| Battery size | 10 | 10 | 312/10 | 13/312 |
| Power consumption | 0.9 mA | 0.9 mA | 1.0 mA | 1.0 mA |
| Features | | | | |
| BestSound Technology | micon | micon | micon | micon |
| e2e wireless 2.0 | - | ○ | ○ | ○ |
| Volume control | - | ○ ⁽⁶⁾ | ○ | ○ |
| Push button | - | ○ | ○ | ○ |
| T-coil | - | - | ○ | ○ |
| Audio input | - | ○ ⁽⁶⁾ | ○ ⁽⁶⁾ | ○ ⁽⁶⁾ |
| Tinnitus noiser | ● | ● | ● | ● |
| Optivent | ● | ● | ● | ● |
| Accessories | | | | |
| Programming adapter | Flex connector | Flex connector | Flex connector | Flex connector |

(6) With miniTek
(11) Max. values


● Standard
○ Optional

Orion.

Flexible, solid, individual essentials.

Orion™ combines smart technology with long-lasting durability for solid hearing aids that provide a reliable listening experience. Its nanocoated housing and IP67-certification combine with improved speech intelligibility and wireless connectivity to make Orion a trustworthy partner – no matter what life has in store for your customers.

Orion. With great comfort and superior connectivity.



For more information, please consult our "Orion" brochure.



Orion

| Type | RIC |
|------------------------------|---|
| Performance level | Essential |
| Description | Orion RIC |
| Picture | |
| Fitting range | |
| Product details | miniReceiver S 45/108 dB miniReceiver M 60/118 dB miniReceiver P 70/124 dB miniReceiver HP 75/130 dB |
| Matrix (2 ccm) | |
| Battery size | 312 |
| Power consumption | S 0.8, M/P 1.0, HP 1.1 mA |
| Features | |
| e2e wireless 2.0 | ● |
| Rechargeable | - |
| Volume control | ○ ¹ |
| Push button | ○ ¹ |
| T-coil | ● |
| Battery door on/off function | ● |
| LifeTube/Standard earhook | - |
| Audio input | - |
| Tinnitus noiser | ● |
| Optivent | - |
| Modifications | |
| Color conversion kit | ● |
| Accessories | |
| Programming cable/adaptor | size 312 |
| Fitting set | miniReceiver 2.0 |

Accessories, see pages 56 – 63.



Flexible color concept

The housings of Orion's S, M and P models consist of the bottom shell in Dark Granite or Beige that can be combined with a selection of exchangeable top shells in a variety of colors.



Colors:

BTE: S/M/P



RIC/BTE: SP



ITE



| BTE | | | | ITE |
|---|---|---|-------------------|---|
| Essential | | | | Essential |
| Orion S | Orion M | Orion P | Orion SP | Orion CIC, ITC, ITE |
| | | | | |
| | | | | |
| LifeTube 45/124 dB Earhook 55/124 dB | LifeTube 53/125 dB Earhook 60/130 dB | LifeTube 65/125 dB Earhook 70/130 dB | Earhook 82/140 dB | CIC 50/113 ¹¹ dB ITC 55/118 ¹¹ dB ITE 60/123 ¹¹ dB |
| 312 0.8 mA | 13 1.0 mA | 13 1.0 – 1.3 mA | 675 2.1 mA | CIC 10, ITC 312/10, ITE 13/312 CIC 0.9 mA, ITC/ITE 1.0 mA |
| ● | ● | ● | ● | ● |
| - | - | - | - | - |
| ○ ⁷ | ○ ⁷ | ● | ● | ○ ¹⁰ |
| ● ⁷ | ○ ⁷ | ● | ● | ○ |
| - | ● | ● | ● | ○ ¹⁰ |
| ● | ● | ● | ● | ● |
| ● | ● | ● | ● | - |
| - | ○ ⁸ | ○ ⁸ | ○ ⁸ | - |
| ● | ● | ● | ● | ● |
| - | - | - | - | ○ |
| ● | ● | ● | - | - |
| size 312 LifeTube | size 13 LifeTube | size 13 LifeTube | CS 44 - | Flex connector - |

- (1) Exchangeable controls – programmable via Connexx
- (7) Programmable via Connexx
- (8) With battery compartment modification
- (10) Not available for CIC
- (11) Max. values

- Standard
- Optional

Sirion.

Simply reliable.

Built for durability and dependability: With a sturdy, IP67-rated housing that protects the technology inside, Sirion™ offers a reliable performance when it comes to speech intelligibility – even in noisy environments. The combination of customizable settings and simple convenience makes Sirion a simple and reliable choice for wearers who don't want to miss out on the action in their lives. **Sirion. With great comfort and versatility.**



For more information, please consult our "Sirion" brochure.

Sirion

| Type | BTE |
|------------------------------|---|
| Performance level | Essential |
| Description | Sirion S |
| Picture | |
| Fitting range | |
| Product details | |
| Matrix (2 ccm) | LifeTube 45/124 dB Earhook 55/124 dB |
| Battery size | 312 |
| Power consumption | 0.8 mA |
| Features | |
| e2e wireless 2.0 | - |
| Rechargeable | - |
| Volume control | ○ ⁷ |
| Push button | ● ⁷ |
| T-coil | - |
| Battery door on/off function | ● |
| LifeTube/Standard earhook | ● |
| Audio input | - |
| Tinnitus noiser | - |
| Optivent | - |
| Modifications | |
| Color conversion kit | ● |
| Accessories | |
| Programming cable/adaptor | size 312 |
| Fitting set | LifeTube |

Microphone

Using the natural ear for directivity, the omnidirectional microphone takes in sound from all directions.

Protective nanocoated housing

Nanocoated membranes for the directional microphones, a closed frame and a rubber gasket prevent moisture, sweat, dust and dirt from entering the hearing aid.


Rotary volume control¹⁰

Vent
For optimum ear ventilation. With Optivent option.



Sirion S
Sirion ITE

Actual size in cm



Flexible color concept

The housings of Sirion's BTE models consist of the bottom shell in Dark Granite or Beige that can be combined with a selection of exchangeable top shells in a variety of colors.



Colors:

BTE



ITE



See details on page 91.

| BTE | | ITE |
|---|---|---|
| Essential | | Essential |
| Sirion M | Sirion P | Sirion CIC, ITC, ITE |
| | | |
| | | |
| LifeTube 53/125 dB Earhook 60/130 dB | LifeTube 65/125 dB Earhook 70/130 dB | CIC 50/113 ¹¹ dB ITC 55/118 ¹¹ dB ITE 60/123 ¹¹ dB |
| 13 | 13 | CIC 10, ITC 312/10, ITE 13/312 |
| 1.0 mA | 1.0 – 1.3 mA | CIC 0.9 mA, ITC/ITE 1.0 mA |
| - | - | - |
| ○ ⁷ | ● | ○ ¹⁰ |
| ○ ⁷ | ● | ○ |
| ● | ● | ○ ¹⁰ |
| ● | ● | ● |
| ● | ● | - |
| ○ ⁸ | ○ ⁸ | - |
| - | - | - |
| - | - | ○ |
| ● | ● | - |
| size 13 | size 13 | Flex connector |
| LifeTube | LifeTube | - |

(7) Programmable via Connexx

(8) With battery compartment modification

(10) Not available for CIC


(11) Max. values

● Standard

○ Optional

Intuis.


Understanding – made easy.



Intuis™ stands for a straightforward, comfortable and reliable line of hearing aids at a good value. Help your customers enjoy a dependable hearing experience with its hassle-free handling, solid and attractive housings and proven hearing system technology in both BTE and custom models. And with the improved audibility and higher sound quality of three Intuis Pro models, you can offer even more of your customers the joy of hearing comfort in any situation. Features like adjustable 4 Channel Adaptive Noise Reduction and Speech Direct reduce listening effort and background noise for more natural, comfortable and high-class hearing in any situation.

Intuis. Easy to wear.

For more information, please consult our "Intuis" brochure.



Intuis Pro

Colors:



| Type | BTE | | |
|---|---|-------------------|-------------------|
| Performance levels | Basic | | |
| Description | Intuis Pro S Dir | Intuis Pro Dir | Intuis Pro SP Dir |
| Picture | | | |
| Fitting range | | | |
| Product details | | | |
| Matrix (2 ccm) | LifeTube 53/123 dB Earhook 55/124 dB | Earhook 60/130 dB | Earhook 80/138 dB |
| Battery size | 13 | 13 | 675 |
| Power consumption | 0.7 mA | 0.8 mA | 1.0 mA |
| Features | | | |
| Frequency channels/ Compression controls | 4/4 | 4/4 | 4/4 |
| Hearing programs | 1 – 4 | 1 – 4 | 1 – 4 |
| e2e wireless 2.0 | - | - | - |
| Rechargeable | - | - | - |
| Volume control | - | ● | ● |
| Push button | ● | ● | ● |
| T-coil | ● | ● | ● |
| Battery door on/off function | ● | ● | ● |
| LifeTube/Standard earhook | ● | ● | ● |
| Small earhook | ● | ● | ● |
| Audio input | - | ● | ● |
| Tinnitus noiser | - | - | - |
| Modifications | | | |
| CROS | ● ⁽¹³⁾ | ● | ● |
| BICROS | ● ⁽¹³⁾ | ● | ● |
| Eyeglass adapter | ○ | ○ | ○ |
| Accessories | | | |
| Programming cable | CS 44 | CS 44 | CS 44 |
| Fitting set | LifeTube | - | - |

(13) With eyeglass modification

● Standard
○ Optional

Colors:

BTE



ITE



Intuis

| Type | BTE | |
|---|---|---|
| Performance level | Basic | |
| Description | Intuis Life | Intuis S Dir |
| Picture | | |
| Fitting range | | |
| Product details | | |
| Matrix (2 ccm) | LifeTube 35/115 dB Earhook 45/117 dB | LifeTube 53/123 dB Earhook 55/124 dB |
| Battery size | 312 | 13 |
| Power consumption | 0.5 mA | 0.7 mA |
| Features | | |
| Frequency channels/ Compression controls | 4/4 | 4/4 |
| Hearing programs | 1 | 4 |
| e2e wireless 2.0 | - | - |
| Rechargeable | - | - |
| Volume control | - | - |
| Push button | - | ● |
| T-coil | - | ● |
| Battery door on/off function | ● | ● |
| LifeTube/Standard earhook | ● | ● |
| Small earhook | - | ● |
| Audio input | - | - |
| Tinnitus noiser | - | - |
| Modifications | | |
| CROS | ● ¹³ | ● ¹³ |
| BICROS | ● ¹³ | ● ¹³ |
| Eyeglass adapter | ○ | ○ |
| Accessories | | |
| Programming cable/adaptor | CS 44 | CS 44 |
| Fitting set | LifeTube | LifeTube |

| BTE | | ITE |
|-------------------|-------------------|---|
| Basic | | Basic |
| Intuis Dir | Intuis SP Dir | Intuis CIC, ITC, ITE |
| | | |
| | | |
| Earhook 60/130 dB | Earhook 80/138 dB | CIC 47/113 ¹¹ dB ITC 50/118 ¹¹ dB ITE 60/123 ¹¹ dB |
| 13 | 675 | CIC 10, ITC 312, ITE 13 |
| 0.8 mA | 1.0 mA | 0.6 mA |
| Features | | |
| 4/4 | 4/4 | 4/4 |
| 4 | 4 | ○ |
| - | - | - |
| - | - | - |
| ● | ● | ○ ¹⁰ |
| ● | ● | ○ |
| ● | ● | ○ ¹⁰ |
| ● | ● | - |
| ● | ● | - |
| ● | ● | - |
| ● | ● | - |
| ● | ● | - |
| ● | ● | - |
| ○ | ○ | - |
| Accessories | | |
| CS 44 | CS 44 | Flex connector |
| - | - | - |

(10) Not available for CIC
(11) Max. values


(13) With eyeglass modification

● Standard
○ Optional

Lotus.

The reliable solution.

Lotus™ combines proven technology with ease of use and good audio quality – at outstanding value. The housing is robust and particularly reliable for every day use. With its broad range of models, Siemens Lotus offers an appropriate and comfortable solution for any wearer. **Lotus. For everything you love to hear.**



Lotus Pro Programmable

| Type | BTE |
|--|-------------------|
| Performance level | Basic |
| Description | Lotus Pro 2 M |
| Picture | |
| Fitting range | |
| Product details | |
| Matrix (2 ccm) | Earhook 60/130 dB |
| Battery size | 13 |
| Power consumption | 0.8 mA |
| Features | |
| Frequency channels/ Compression controls | 2/1 |
| Hearing programs | 3+A |
| e2e wireless 2.0 | - |
| Rechargeable | - |
| Volume control | ● |
| Push button | ● |
| T-coil | ● |
| Battery on/off function | ● |
| LifeTube/Standard earhook | ● |
| Small earhook | ● |
| Audio input | ● |
| Tinnitus noiser | - |
| Modifications | |
| CROS | ● |
| BICROS | ● |
| Eyeglass adapter | ● |
| Accessories | |
| Programming cable/adaptor | CS 44 |
| Fitting set | - |



For more information, please consult our "Lotus" brochure.

Colors:

BTE



ITE



| BTE | | | | ITE |
|-------------------|---|-------------------|-------------------|---|
| Basic | | | | Basic |
| Lotus Pro 2 SP | Lotus Pro M | Lotus Pro P | Lotus Pro SP | Lotus Pro CIC, ITC, ITE |
| | | | | |
| | | | | |
| Earhook 80/138 dB | LifeTube 47/126 dB Earhook 60/130 dB | Earhook 70/133 dB | Earhook 80/138 dB | CIC 47/113 ¹¹ dB ITC 50/118 ¹¹ dB ITE 60/123 ¹¹ dB |
| 675 | 13 | 13 | 675 | CIC 10, ITC 312, ITE 13 |
| 1.6 mA | 0.8 – 1.0 mA | 0.7 mA | 2.0 mA | 0.6 mA |
| 2/1 | 2/1 | 2/1 | 2/1 | 2/1 |
| 3+A | 3+A | 2+A | 2+A | ○ |
| - | - | - | - | - |
| - | - | - | - | - |
| ● | ● | ● | ● | ○ ¹⁰ |
| ● | ● | ● | ● | ○ ¹⁰ |
| ● | ● | ● | ● | ○ ¹⁰ |
| ● | ● | ● | ● | - |
| ● | ● | ● | ● | - |
| ● | - | ● | ● | - |
| ● | ● | ● | ● | - |
| - | - | - | - | - |
| ● | ● | ● | ● | - |
| ● | ● | ● | ● | - |
| ● | ○ | ● | ● | - |
| CS 44 | CS 44 | CS 44 | CS 44 | Flex connector |
| - | LifeTube | - | - | - |

(10) Not available for CIC
(11) Max. values

● Standard
○ Optional

Lotus Trimmer

Colors:

BTE



ITE



| Type | BTE | | |
|------------------------------|-------------------|-------------------|-------------------|
| Performance level | Basic | | |
| Description | Lotus 23 M | Lotus 23 P | Lotus 23 SP |
| Picture | | | |
| Fitting range | | | |
| Product details | | | |
| Matrix (2 ccm) | Earhook 60/130 dB | Earhook 70/132 dB | Earhook 80/138 dB |
| Battery size | 13 | 13 | 675 |
| Power consumption | 1.0 mA | 0.9 mA | 0.9 mA |
| Features | | | |
| Frequency channels | 2 | 2 | 2 |
| Trimmers | NH/MPO/Gain | NH/MPO/Gain | NH/MPO/Gain |
| Hearing programs | 2 | 2 | 2 |
| e2e wireless 2.0 | - | - | - |
| Rechargeable | - | - | - |
| Volume control | ● | ● | ● |
| Push button | ● | ● | ● |
| T-coil | ● | ● | ● |
| Battery door on/off function | ● | ● | ● |
| Standard earhook | ● | ● | ● |
| Small earhook | - | ● | ● |
| Audio input | ● | ● | ● |
| Tinnitus noiser | - | - | - |
| Modifications | | | |
| CROS | ● | ● | ● |
| BICROS | ● | ● | ● |
| Eyeglass adapter | - | - | - |
| Accessories | | | |
| Programming cable/adaptor | CS 44 | CS 44 | CS 44 |

(10) Not available for CIC
(11) Max. values

(14) With screw set volume control for CIC

| ITE | | | |
|------------------------------|---------------------------|---------------------------|---------------------------|
| Basic | | | |
| Description | Lotus 23 CIC | Lotus 23 ITC | Lotus 23 ITE |
| Picture | | | |
| Fitting range | | | |
| Product details | | | |
| Matrix (2 ccm) | 40/113 ⁽¹¹⁾ dB | 50/118 ⁽¹¹⁾ dB | 60/123 ⁽¹¹⁾ dB |
| Battery size | 10 | 312 | 13 |
| Power consumption | 0.7 mA | 0.5 mA | 0.5 mA |
| Features | | | |
| Frequency channels | 2 | 2 | 2 |
| Trimmers | NH/MPO/NR | NH/MPO/NR | NH/MPO/NR |
| Hearing programs | 1 | 1 | 1 |
| e2e wireless 2.0 | - | - | - |
| Rechargeable | - | - | - |
| Volume control | Manual VC ⁽¹⁴⁾ | - | - |
| Push button | - | - | - |
| T-coil | ○ | ○ | ○ |
| Battery door on/off function | - | - | - |
| Standard earhook | - | - | - |
| Small earhook | - | - | - |
| Audio input | - | - | - |
| Tinnitus noiser | - | - | - |
| Modifications | | | |
| CROS | - | - | - |
| BICROS | - | - | - |
| Eyeglass adapter | - | - | - |
| Accessories | | | |
| Programming cable/adaptor | Flex connector | Flex connector | Flex connector |

● Standard
○ Optional

SecureEar.

Help your customers experience more by sharpening their senses.



SecureEar™ is an active digital hearing protection that prevents the wearer from excessive noise and loud sounds. It ideally combines hearing protection with the ability to communicate effectively in all environments by using many of the technologies available in the portfolio of Siemens hearing aids. Customization ensures a perfect fit. Whatever their environment, e.g. a day at the shooting range or working with heavy machinery, your customers will be protected, yet still able to hear what's going on around them. With SecureEar, the wearer can also receive cell phone or CB calls at any time thanks to the integrated T-coil.

SecureEar with T-coil: Protects hearing and also allows the user to communicate without any constraints.



Protection features

- Hearing protection with 85 dB maximum output sound
- Tested according to EN 352-2:2002 and EN 352-7:2002

Communication features

- 4 channel sound processing for enhanced hearing and communication
- Noise reduction
- Volume control with 24 dB range
- Maximal amplification 30 dB
- Optional T-coil

Sound attenuation

| Frequency (Hz) | Sound attenuation (dB) | Standard deviation (dB) | APV (dB) |
|----------------|------------------------|-------------------------|----------|
| 63 | 29.9 | 4.9 | 25.0 |
| 125 | 34.0 | 6.2 | 27.8 |
| 250 | 33.0 | 4.6 | 28.4 |
| 500 | 33.7 | 4.8 | 28.9 |
| 1000 | 37.1 | 3.5 | 33.6 |
| 2000 | 40.1 | 2.2 | 37.9 |
| 4000 | 42.6 | 4.4 | 38.2 |
| 8000 | 43.1 | 4.9 | 38.2 |

SNR = 36 dB, H = 37 dB, M = 33 dB, L = 30 dB

APV – Assumed Protection Value, SNR – Single Number Rating attenuation, H – High-frequency attenuation, M – Medium-frequency attenuation, L – Low-frequency attenuation

SecureEar

Colors:

Faceplate

 Beige Tan

Shell

 Beige Tan Brown Blue Red Orange

| | |
|--|--------------------|
| Type | ITE |
| Description | Hearing protection |
| Picture | |
| Output sound pressure level | |
| Active sound pressure | |
| Product details | |
| Matrix (2 ccm) | 85 dB |
| Battery size | 312 |
| Battery life | ~ 230 h |
| Features | |
| Volume control | ● |
| T-coil | ○ |
| Housing/Shell | |
| A perfect fit. The shell is made of robust UV material. | |
| Certification | |
| EN 352-2:2002 EN 352-7:2002 | |

● Standard
○ Optional

Aloe Vera gel

This gel is a skin care gel. The gel uses Aloe Vera as its active ingredient to prevent irritation, condition the skin, and make it easier to insert SecureEar. Use of the gel is recommended before and after wearing SecureEar to ensure continued wearing comfort.



Direct audio input and child-lock battery door.

The added value of individual options: Special input and safety features.

Meet the needs of each customer with flexible options for individual requirements. Many of our world-class hearing aids are available with options that give you precisely this flexibility. Features like optional audio shoes give patients access to direct audio input, which can make a world of difference when it comes to improving speech understanding in noisy situations. And the optional child-lock battery door can protect young wearers and patients who require assistance.

Cutting through the noise: Direct audio input.

Give wearers the ability to hear speakers clearly, even in noisy environments. Optional direct audio input enables understanding by connecting to FM systems. It is particularly useful for children so that they can hear their teachers over the noise of the classroom, or for wearers attending meetings so that they can focus on the speaker.

Keeping wearers safe: Child-lock battery doors.

Give parents and caretakers peace of mind by offering hearing instruments that include safety options like the child-lock battery door. This feature keeps batteries securely in place, reducing the risk of damage or choking on small parts.



The standard battery compartment of Motion SA binax models can be exchanged for a solution with a built-in audio shoe.



Motion M and P micon models are equipped with a child-lock battery door, so children don't unintentionally open it.

| Products | Colors and Item No. | |
|---|--|-------------------------------------|
| Carat A binax | | |
| Battery door with built-in audio shoe | Beige: 10822542 / Dark Granite: 10823652 | |
| Motion SA binax | | |
| Battery door with built-in audio shoe | Beige: 10822542 / Dark Granite: 10823652 | |
| Motion M micon and Motion P micon | M | P |
| Pediatric kit including child-safe battery door and accessories | Beige / Silver: 10820655 | Beige / Silver: 10820636 |
| Battery door with child-lock | Beige: 10668306 Silver: 10821034 | Beige: 10676204 Silver: 10820538 |
| Audio shoe transparent | 10820642 | 10820642 |
| Nitro micon and Orion SP | | |
| Pediatric kit including child-safe battery door and accessories | Beige: 10820940 / Granite: 10821160 / Grey: 10821157 / Brown: 10821159 Black: 10821158 / Dark Granite: 10821156 | |
| Pediatric service battery door | Beige: 10822106 / Dark Granite: 10822107 / Grey: 10822108 Black: 10822109 / Brown: 10822110 / Granite: 10822111 | |
| Orion M/P and Sirion M/P | | |
| Battery door with built-in audio shoe | Beige: 10673021 / Dark Granite: 10674609 | |
| Tamperproof battery door | Beige: 10820326 / Dark Granite: 10820333 | |

Carat A binax

Premium binaural technology, easy handling and compact size, all with optional direct audio input.

Other features include:

- IP67-rated for resistance to moisture, sweat, dust and dirt
- Standard battery compartment can be exchanged for a DAI-ready version with built-in audio shoe
- Exceptional wireless connectivity



Carat A binax

Performance levels

- 7bx
- 5bx
- 3bx **New!**

For more details, see page 22.

Motion SA binax **New!**

Premium binaural technology, compact design, and fully featured, all with optional direct audio input.

Other features include:

- IP67-rated for resistance to moisture, sweat, dust and dirt
- Standard battery compartment can be exchanged for a DAI-ready version with built-in audio shoe
- Sleek housing: Motion SA binax is almost imperceptible when worn



Motion SA binax

Performance levels

- 7bx
- 5bx
- 3bx

For more details, see page 24.

Motion M and P micon

BestSound Technology, optimized design, and many safety options provide outstanding functionality.

Other features include:

- IP67-rated for resistance to moisture, sweat, dust and dirt
- DAI-ready with optional audio shoe
- Programmable key lock and child-lock battery door as options



Motion P micon

Motion M micon

Performance levels

- 7mi
- 5mi
- 3mi

For more details, see page 32.

Nitro micon

Maximum amplification in all modes makes this super power hearing instrument shine.

Other features include:

- IP67-rated for resistance to moisture, sweat, dust and dirt
- Long-lasting 675 size standard battery
- Optional child-lock battery door with audio input option



Nitro micon

Performance levels

- 7mi
- 3mi

For more details, see page 36.

Orion/Sirion M and P models

For extra security, the standard battery door of Orion and Sirion M and P models can be replaced with a tamperproof battery door so the battery compartment cannot be opened unintentionally. The standard battery door can also be exchanged with a battery door with design-integrated audio shoe.

Other features include:

- IP67-rated for resistance to moisture, sweat, dust and dirt
- Robust housing



Orion P
Sirion P

Orion M
Sirion M

Performance levels

- Essential

For more details, see page 40–43.



Accessories and hearing aid compatibility.


This overview shows you which accessories work with Siemens binax and micon hearing instruments. Depending on the solution your customers choose, they can even enhance their hearing experience and convenience with these well-matching audio streamers, audio streamer accessories, remote controls, and charging stations.

Useful Accessories.

Audio streaming

easyTek and **miniTek™** audio streamers work seamlessly with our hearing aids. These modern accessories wirelessly connect the hearing instruments to many different devices. Sound from Bluetooth phones, mp3 players, televisions and other audio sources is streamed directly to the hearing aids. Plus, they also act as remote controls for discreet hearing program and volume adjustments. The smartphones can be used to conveniently operate the **easyTek** via the **easyTek App***, and the **miniTek** via the **miniTek Remote App****. The apps can be downloaded on the Google Play Store or Apple App Store at no cost.

The **VoiceLink™** companion microphone connects to easyTek and miniTek and effectively transmits speech when wearers are attending conferences or presentations.



Remote controls

The **touchControl App*** for binax hearing aids lets your clients discreetly and conveniently change their hearing programs, volume, bass and treble via their smartphones – without any extra hardware. The sleek **easyPocket™** with large buttons and an easy-to-read display gives even users with dexterity issues easy access to their hearing instruments.



Charging station

Siemens **eCharger** recharges and dehumidifies hearing aids overnight, so your clients don't need to worry about changing tiny batteries.




| Hearing aids | Audio streamers | Audio streamer accessories | Remote controls | Charging station |
|---------------------|-----------------|--|-----------------------------|------------------|
| | easyTek miniTek | easyTek App miniTek Remote App Transmitter | touchControl App easyPocket | eCharger |
| Ace binax | - | - | • | - |
| Pure binax | • | • | • | • |
| Carat binax | • | • | • | • |
| Carat A binax | • | • | • | - |
| Motion SX, PX binax | • | • | • | • |
| Motion SA binax | • | • | • | - |
| Insio binax | • | • | • | - |
| Pure micon | - | • | - | • |
| Siemens Life micon | - | • | - | • |
| Motion SX, PX micon | - | • | - | • |
| Motion M, P micon | - | • | - | • |
| Aquaris micon | - | • | - | • |
| Nitro micon | - | • | - | • |
| Insio micon | - | • | - | • |
| Orion | - | - | - | • |

* Compatible with iOS and Android.
** Compatible with Android.


All possibilities in your clients' hands: easyTek.

Satisfy even your most tech-savvy customers: Worn around the neck, the lightweight easyTek* audio streamer and remote control quickly and reliably connects almost all binax hearing aids to Bluetooth-enabled devices. Ready right out of the box, it feeds sound from TVs, mp3 players, phones, the VoiceLink companion microphone, or other audio sources directly to wearers' instruments. On our website, you can watch useful handling videos and find out more about this high-tech device.

* Compatible with Android and iOS.



Audio streaming:




easyTek connects Pure binax, Carat binax, Carat A binax, Motion binax and Inso binax to iPhones, iPads, Android devices and to any Bluetooth-enabled phones.



easyTek

easyTek – easy pairing.

Before easyTek can be used, it needs to be paired with several devices. And easyTek's name says it all – pairing is as easy as 1, 2, 3.




Hearing aids



Step 1:

- Put easyTek around your customer's neck.
- Turn easyTek on by pressing **1** until the green LED **5** turns on.



Step 2:

- Put the hearing instruments into the ears and turn them on.
- The green LED **5** flashes briefly and then slowly fades.
- Now the hearing aids are paired.

Bluetooth-enabled devices



Step 1:

- Turn easyTek on by pressing and holding **1** for several seconds.
- Press **1** and **2** simultaneously until the blue LED **4** flashes.



Step 2:

- Turn on the device's Bluetooth function.
- Search for easyTek and select it.
- The blue LED **4** flashes briefly and then slowly fades.

Transmitter



Step 1:

- Turn easyTek on by pressing and holding **1** for several seconds.
- Press **1** and **2** and **3** simultaneously until the blue LED **4** and the green LED **5** flash.



Step 2:

- Turn the transmitter on within 1 m of easyTek and wait for the pairing to complete. This takes approximately 2 minutes.
- The blue LED **4** will stop flashing and remain on.

Further functionality: The easyTek App.

Even exceeding easyTek in functionality, the easyTek App* is an ultra-discreet, easy way to control binax hearing aids. Wearers simply look like they are checking their smartphones, yet can quickly control hearing programs, Spatial Configurator, volume, tinnitus settings, SoundBalance™ – even while streaming music.

The easyTek App can be downloaded for free on the Google Play Store or Apple App Store.



Full control at your clients' fingertips: touchControl App.

The touchControl App lets your clients discreetly and conveniently change their hearing programs, volume, bass and treble via their Android or iOS smartphones. With an operating range of 1 meter, it sends an ultra-high frequency signal directly to the hearing aids – without needing any extra hardware. The touchControl App works with all binax hearing aids and is free to download on the Google Play Store or Apple App Store.

Unprecedented discretion: Even the smallest binax hearing aids are remotely controllable via the touchControl App.



Insio binax IIC Ace binax

Display current program

SoundBalance for bass and treble adjustments

Volume up/down

Mute

Further options

Easy control of binax hearing aid functions via the user's smartphone.

Main screen

Further options

Available on the iPhone
App Store

ANDROID APP ON
Google play

binax apps.

Added value for your clients – and for you.

Delighting wearers and providing remarkable value is essential in today's hearing care market. Help your clients download the easyTek App and touchControl App on the Google Play Store or Apple App Store. It's so simple – and for free!



App download.

- Make sure your customer has an Android or iOS smartphone and an account on the Google Play Store or Apple App Store.
- Together with your customer, log into their account on the Google Play Store or Apple App Store with their Android or iOS smartphone. Find the easyTek App or touchControl App.
- Follow the instructions until the app is installed on your customer's Android or iOS smartphone.
- Show your customer how easy it is to control their hearing instruments via the respective app.

You can also find a link to the app stores on the following webpage

www.bestsound-technology.com/apps






Audio Streamers.

| Products | Type | Description | Item No. |
|---|------------------------|--|---|
| easyTek | | | |
|  | easyTek | easyTek audio streamer and remote control incl. neckloop antenna. Transmitters and Voicelink need to be purchased separately. Size: 49x12 mm | 10824190* 10931133** |
|  | easyTek App | The easyTek App allows Android and iOS users to control their easyTek via their smartphone. The interface has been carefully designed so it is easy to understand. | Can be downloaded on the Google Play Store or Apple App Store at no cost. |
| miniTek | | | |
|  | miniTek set | miniTek, transmitter, lanyard, power supply set 2 x, accessory set Size: 61x45x15,5 mm | 10603493 |
|  | 1. miniTek pouch black | | 10604181 |
| | miniTek pouch grey | | 10604619 |
| | 2. Car adapter | For Tek/miniTek and eCharger | 10601544 |
| | 3. Transmitter | | 10603245 |
| | 4. Belt clip | | 10604180 |
| | 5. Lanyard | | 10604232 |
| | Power supply set | | 10417270 |
| | Programming cable | USB cable for Connexx | 10417279 |
| | Cable set | | 10420199 |
| | Adapter set | | 10420201 |
| Slip case | | 10603494 | |
|  | miniTek Remote App | The miniTek Remote App allows Android users to control their miniTek via their smartphone. And because the app's interface has the same look as the miniTek, it is easy to understand. | Can be downloaded on the Google Play Store at no cost. |
| Voicelink | | | |
|  | Voicelink set | Transmitter, lapel microphone, power supply set, audio cable set, adapter set | 10673506 |
| | Microphone service set | Lapel microphone, accessory set | 10673511 |

* Compatible with Pure binax, Motion binax, Carat & Carat A binax.
** Compatible with Insio binax and other binax products featuring e2e wireless 3.0.

Remote Controls and eCharger.

| Products | Type | Description | Item No. |
|--|---|---|---|
| touchControl App | | | |
|  | touchControl App | The touchControl App allows Android and iOS users to control their hearing instrument programs, volume, bass and treble via their smartphone. | Can be downloaded on the Google Play Store or Apple App Store at no cost. |
| easyPocket | | | |
|  | easyPocket | User-friendly, ergonomic remote control with a modern design and built-in display. For wireless enabled hearing aids 7bx, 5bx, 3bx, 7mi, 5mi, 3mi, Orion. Size: 92,5x42,5x17,5 mm | 10675287 |
| eCharger | | | |
|  | eCharger* with electronic drying function | eCharger without inlay | 10600230 |
| | Inlay 312-B | Inlay for Pure XCEL**, X01, X00 | 10423516 |
| | Inlay 13-C | Inlay for Motion X00 M, M VC, P Motion X01D, D M VC, D P | 10423515 |
| | Inlay 13-E | Inlay for Motion X00 S | 10423076 |
| | Inlay 13-F | Inlay for Motion XCEL P**, X01 P | 10600074 |
| | Inlay 13-G | Inlay for Motion SX micon, XCEL SX**, X01 SX | 10600297 |
| | Inlay 13-H | Inlay for Pure XCEL Carat, Pure Carat | 10603039 |
| | Inlay 13-K | Inlay for Motion PX binax, Motion PX micon | 10823023 |
| | Inlay 13-M | Inlay for Carat binax | 10824799 |
| | Inlay 13-N New! | Inlay for Motion SX binax | 10827814 |
| | Inlay 312-L | Inlay for Pure binax, Pure micon | 10820446 |
| | Rechargeable battery 312 | For size 312 rechargeable hearing aids | 10423863 |
| | Rechargeable battery 13 | For size 13 rechargeable hearing aids | 10423862 |
| | Car adapter | For eCharger, miniTek and Tek | 10601544 |
| | KuanTech power supply | With 4 power supply adapters | 10666693 |

* Rechargeable battery not included. ** Not for performance level 101 XCEL.

miniReceiver 2.0

More robust, more flexible.

The miniReceiver 2.0 is a smart RIC receiver that uses a smaller version of SIM technology. It securely and reliably stores and transmits data between the receiver and the hearing aid chip. The SIM card can also detect an incorrect receiver and automatically reduce gain to protect the wearer. And the design enables more ergonomic placement of the receiver for better microphone placement. This improves directionality, reduces feedback sensitivity, and expands the fitting range – for better sound quality, enhanced audiological performance and improved usability.

miniReceiver 2.0 – ready to perform each and every day.

Robust connector
Protects against dust, dirt, and moisture.

Reliable click interface
Proven connection to Click Domes™ and Click Molds™.

Effective 2-step wax protection
In combination with Click Domes and Click Molds.

SIM Card Smart technology


- Auto identification safety feature
- Reliable and secure logical connection
- Enables better sound quality

- 1 Front microphone moved to top
- 2 SIM card
- 3 Plug and jack moved down

Color indicator
Blue left side and red right side immediately identifiable.

HF4 Pro filter New!
The new high-tech filter gauze reliably protects the hearing instrument by repelling dirt, cerumen, oil and water – for an extra-long usability.

miniReceiver 2.0 – exchanging is as easy as 1, 2, 3




1 Open the locking bar. 2 Remove/Insert the miniReceiver 2.0. 3 Close the locking bar.

The miniReceiver 2.0 is easily exchangeable. To replace the miniReceiver 2.0, use the exchange tool to open and close the locking bar for reliable securing and no risk of dislodgement.

To see an animation of how to exchange the miniReceiver 2.0 visit the website for professionals.










miniReceiver Fitting.

For every generation of BestSound Technology, we offer a wide range of miniReceivers for every requirement.




| Products | Description | Item No. | | | | |
|--|--|--|---|------------------|----------|----------|
| miniReceiver Fitting Set | miniReceiver 2.0 Fitting set miniReceiver S units (1 set of sizes 1–3), miniReceiver M/P units (1 set of size 2), Click Dome packs (1 of each size), Click Mold service set, filter HF4 left and right (15 filters of each), measuring card for receiver length, measuring card for Click Mold size | 10827741 | | | | |
| miniReceiver Type | Size | Units/box | miniReceiver | miniReceiver 2.0 | | |
| Type S miniReceiver • Pure XCEL • Pure XCEL Carat • Pure • Pure Carat • Pure SE | miniReceiver 2.0 • Ace binax • Pure binax • Carat binax • Ace micon • Pure micon • Orion RIC | 0L | 1 | 10425988 | 10669907 | |
| | | 0R | 1 | 10425989 | 10669908 | |
| | | 1L | 1 | 10425990 | 10669909 | |
| | | 1R | 1 | 10425991 | 10669910 | |
| | | 2L | 1 | 10425992 | 10669911 | |
| | | 2R | 1 | 10425993 | 10669912 | |
| | | 3L | 1 | 10425994 | 10669913 | |
| | | 3R | 1 | 10425995 | 10669914 | |
| | Type M miniReceiver • Pure XCEL • Pure XCEL Carat • Pure • Pure Carat | miniReceiver 2.0 • Ace binax • Pure binax • Carat binax • Ace micon • Pure micon • Orion RIC | 0L | 1 | 10671603 | 10669915 |
| | | | 0R | 1 | 10671604 | 10669916 |
| | | | 1L | 1 | 10671605 | 10669917 |
| | | | 1R | 1 | 10671606 | 10669918 |
| | | 2L | 1 | 10671607 | 10669919 | |
| | | 2R | 1 | 10671608 | 10669920 | |
| | | 3L | 1 | 10671616 | 10669921 | |
| | | 3R | 1 | 10671619 | 10669922 | |
| Type P miniReceiver • Pure XCEL • Pure XCEL Carat • Pure • Pure Carat | | miniReceiver 2.0 • Ace binax • Pure binax • Carat binax • Ace micon • Pure micon • Orion RIC | 0L | 1 | 10426004 | 10669923 |
| | | | 0R | 1 | 10426005 | 10669924 |
| | | | 1L | 1 | 10426006 | 10669925 |
| | | | 1R | 1 | 10426007 | 10669926 |
| | | 2L | 1 | 10426008 | 10669927 | |
| | | 2R | 1 | 10426009 | 10669928 | |
| | | 3L | 1 | 10426010 | 10669929 | |
| | | 3R | 1 | 10426011 | 10669930 | |
| | Type HP miniReceiver • Pure XCEL • Pure XCEL Carat • Pure Carat | miniReceiver 2.0 • Ace binax • Pure binax • Carat binax • Ace micon • Pure micon • Orion RIC | Please order the miniReceiver size HP with a custom CIC shell in the same way as a custom hearing aid using the order form. | | | |










Click Domes and Click Molds.

| Products | Description | Item No. |
|---|---|---|
| Click Dome Types | | |
|  | The open dome is for patients who require maximum comfort and no occlusion. | 4 mm 6 10426024 |
| | | 8 mm 6 10426022 |
| | | 10 mm 6 10426023 |
|  | The semi-open dome combines the benefits of an open and a closed dome for patients who desire the comfort of an open fitting but need more power. | 8-12 mm 6 10669971 |
|  | The closed dome is for patients who require more loudness than the open dome can offer. | 4 mm 6 10426025 |
| | | 8 mm 6 10426020 |
| | | 10 mm 6 10426021 |
|  | The double dome is for patients who require broadband amplification and more retention. Can also be used with a miniReceiver size P. | 8/10 mm 6 10426026 10/12 mm 6 10426027 |
| Click Mold Types | | |
|  | Click Molds are an alternative to fitting with Click Domes. These molds are made by a certified partner. If you have any questions, please do not hesitate to contact us. | Short mold |
|  | | Open mold, short |
|  | | Long mold |
| | | Open mold, long |
| Click Mold User Set | | |
| | Short Click Mold disassemble tool, waxpick and magnet, receiver mold tube set small | 10297177 |
| | Long Click Mold disassemble tool, receiver long mold disassemble set | 10420833 |
| Click Mold Starter Kit | | |
|  | Click Mold Starter Kit Package contents: Click Mold adapters (red and blue), Click Mold service set, Click Mold disassembly tool, Click Mold assembly tool, long mold drilling tool, molding cutters, removal strings, miniReceiver S/M/P dummies, Click Mold measuring card, HF4 filters left and right, miniReceiver wax guard set, work instructions on CD | 10426787 |
| | Click Mold Starter Kit for miniReceiver 2.0 | 10670820 |
| | Click Mold measuring gauge | 10598784 |
|  | Click Mold measuring gauge for miniReceiver 2.0 | 10670821 |

The new HF4 Pro filter.

Better protection for a longer usability: Thanks to its innovative high-tech gauze, the new HF4 Pro filter reliably repels dirt, cerumen, water and oil from hearing aids. Compared to the traditional filter, the interval until the new HF4 Pro filter is clogged and needs to be changed is doubled – and receiver related repairs are significantly reduced.



| Products | Description | Item No. |
|---|--|----------|
| Accessories | | |
|  | 1. HF4 Pro filters right (red) and left (blue) (20 pcs) New! | 10827324 |
| | 2. HF4 filters left (blue) (20x15 pcs) | 10600985 |
| | 3. HF4 filters right (red) (20x15 pcs) | 10600984 |
|  | miniReceiver wax guard (10 pcs) | 10597971 |
|  | 1. miniReceiver 2.0 removal tool | 10674370 |
|  | 2. Click Dome removal tool | 10604944 |
|  | Click Mold service set | 10597930 |
|  | Template for receiver length | 10174571 |
|  | 1. miniReceiver 2.0 measuring gauge For receiver length & impression measurement for P/HP receivers | 10673959 |
|  | 2. miniReceiver HP measuring gauge For receiver length & impression measurement for HP receivers | 10670257 |
|  | 3. Nitro micon measuring gauge The measuring gauge supports the fitting process by letting you easily find out about the right size for the wearer's custom hearing aids. | 10823745 |

The new ThinTube.

New!

Slim shape, secure fit, superb performance.

The new ThinTube is an extremely individual and unobtrusive solution for almost any ear. Its new design is slimmer than previous models and offers a higher cosmetic appeal by using ultra-transparent materials – all while maintaining the same supreme performance and stability. The sleek ThinTube is robust and comfortable, and provides wearers with improved directivity and the amplification they need when rising to the challenges of everyday life.
ThinTube – for a secure fit and optimum sound transmission.



ThinTube
 Ultra-slim and ultra-transparent, yet stable and robust. Equipped with honey tip and concha lock, the ThinTube works together with all proven LifeTips.

Transparent connector
 Firmly attached to the housing. The interface with a screw thread ensures a secure fit between the hearing aid and the ThinTube or the earhook.

LifeTip
 Fits and works perfectly with the new ThinTube.

Designed for a high cosmetic appeal:
 A slimmer design and a more transparent material make the ThinTube more discreet and attractive than previous solutions.




ThinTube

Easy to exchange.

The new ThinTube has a threaded screw connector that is easy to attach and remove from Motion SX and SA binax' new housing. This ensures both a secure fit and a simple exchange process in just a few easy steps.

Thanks to the interface with a screw thread, the ThinTube can be quickly and easily screwed on or off.




ThinTube: Exchange process.



1 Move the ThinTube close to the interface on the hearing aid. Make sure the ThinTube is straight on the screw thread.




2 Rotate the ThinTube clockwise until it sits fully tightened in the screw thread.



3 Rotate the ThinTube until it sits firmly in the interface.




4 Please check if the ThinTube is aligned to the hearing aid body to ensure a perfect fit and optimum sound transmission.



5 To remove the ThinTube, rotate it counterclockwise until it is fully detached from the housing.

LifeTube, S-LifeTube, ThinTube.

| Products | LifeTube |
|----------|---|
| | Compatible hearing aids: Aquaris micon, Intuis, Lotus Pro M, Motion SX micon, Motion XCEL, Motion X01, Orion S/M/P, Siemens Life micon and Sirion S/M/P |

| LifeTube and ThinTube Fitting Sets | Content | Item No. |
|--|--|----------|
|  New! | LifeTubes in all lengths, LifeTips in all sizes, LifeTube accessory box incl. measuring gauge, LifeTube shaper, cleaning wires | 10296968 |

| LifeTube and ThinTube Refill Boxes | Type | Units/box | Item No. |
|--|------|-----------|----------|
|  New! | 1L | 10 | 10296944 |
| | 1R | 10 | 10296945 |
| | 2L | 10 | 10296946 |
| | 2R | 10 | 10296947 |
| | 3L | 10 | 10296948 |
| | 3R | 10 | 10296949 |
| | 4L | 10 | 10296950 |
| | 4R | 10 | 10296951 |
| | 5L | 10 | 10296952 |
| | 5R | 10 | 10296953 |
| | 6L | 10 | 10296954 |
| | 6R | 10 | 10296955 |

| S-LifeTube |
|---|
| Compatible hearing aids: Motion M/P/PX micon, Motion PX binax |

| Content | Item No. |
|--|----------|
| S-LifeTubes in all lengths, LifeTips in all sizes, LifeTube accessory box incl. measuring gauge, LifeTube shaper, cleaning wires | 10296992 |

| Type | Units/box | Item No. |
|------|-----------|----------|
| 1L | 10 | 10297010 |
| 1R | 10 | 10297011 |
| 2L | 10 | 10297012 |
| 2R | 10 | 10297013 |
| 3L | 10 | 10297014 |
| 3R | 10 | 10297015 |
| 4L | 10 | 10297016 |
| 4R | 10 | 10297017 |
| - | - | - |
| - | - | - |
| - | - | - |
| - | - | - |

| ThinTube New! |
|---|
| Compatible hearing aids: Motion SA/SX binax |

| Content | Item No. |
|---|----------|
| ThinTubes in all lengths, measuring gauge, cleaning wires | 10828299 |





| Type | Units/box | Item No. |
|----------------|-----------|----------|
| 1L | 5 | 10828291 |
| 1R | 5 | 10828292 |
| 2L | 5 | 10828293 |
| 2R | 5 | 10828294 |
| 3L | 5 | 10828295 |
| 3R | 5 | 10828296 |
| 4L | 5 | 10828297 |
| 4R | 5 | 10828298 |
| 5 side neutral | 5 | 10828971 |
| - | - | - |
| - | - | - |
| - | - | - |

LifeTips.


| LifeTip Refill Boxes | Type | Units/box | Item No. |
|--|---------------------|-----------|----------|
|  | open 4 mm | 10 | 10296987 |
| | open 6 mm | 10 | 10296988 |
| | open 8 mm | 10 | 10296989 |
| | open 10 mm | 10 | 10296990 |
| | closed 8 mm | 10 | 10296986 |
| | double 8/10 mm | 10 | 10296984 |
| | double 10/12 mm | 10 | 10296985 |
| | semi-open 8 – 12 mm | 10 | 10673779 |

| Additional Accessories | Length | Units/box | Item No. |
|------------------------|--------|-----------|----------|
| LifeTube shaper | 1 | 5 | 10054903 |
| | 2 | 5 | 10054904 |
| | 3 | 5 | 10054905 |
| | 4 | 5 | 10148704 |
| Life cleaning wires | - | 30 | 10054970 |
| Measuring gauge | - | 1 | 10054711 |

Earmolds.

| Earmolds for classic fitting | | | |
|---|---|---|---|
|  |  |  |  |
| Skeleton | Half shell, caved | Full shell, caved | Tube holder |

We offer a comprehensive range of acrylic earmolds for all BTE models. Selected styles are also available as soft earmolds: Skeleton, Canal, Full shell, Half shell.





More comprehensive consultation

Updated information about the full range of Siemens hearing solutions and their benefits helps Hearing Care Professionals give patients clear understanding of all their options.

Detailed information on binaural hearing

Appealing animations and informative explanations let you vividly demonstrate the differences between normal hearing, hearing loss and binaural hearing with hearing aids to your customers.



Performance in detail

Counseling Suite contains comprehensive information on both performance levels and features – with upselling arguments by pointing out the advantages of higher performance levels.

The Counseling Suite.

Competent and flexible consultation – at your fingertips.


Siemens Counseling Suite is a practical tool that facilitates well-structured, customized consultations. Suitable for use with Android and iOS tablets, Windows or Mac PCs, you can use it whenever and wherever you need it.

Intuitive and informative



Hearing Care Professionals can easily tailor the information, animations, and descriptions that they present to each customer with the Counseling Suite's intuitive new design and display options. Its contents give clients a clearer understanding of their hearing loss or tinnitus, its effects, and the benefits of various solutions. For example, the interactive "question and answer" format easily facilitates this process. And the proven hearing loss simulator gives clients and their loved ones both visual and audible insight into their hearing loss.

Individuality

The Counseling Suite lets customers decide which listening situations are most important – and helps them understand the benefits of higher performance levels in those situations.



The Counseling Suite can be used on PCs and tablets – enhancing the flexibility of your customer service. Just download it on the Google Play Store or Apple App Store free of charge!

| Products | Description | Item No. |
|---|---|---|
|  | <p>Counseling Suite The comprehensive advisory tool with practical functions.</p> <ul style="list-style-type: none"> • Helpful support during the consultation process with customers – before, during and after the fitting • Assistance in demonstrating and familiarizing customers with the topics of "hearing" and "hearing instruments," e.g. the hearing loss simulator, selection of the right performance level, maintenance tips • With content on tinnitus, binaural hearing and improved upselling support • Intuitive design for easy use | <p>Available for download at www.bestsound-technology.com/counselingsuite at no cost.</p> |
|  | <p>Counseling Suite App</p> | <p>Available on the Google Play Store or Apple App Store at no cost.</p> |

Connex 7.

More intuitive user interface,
all flexibility.



Connex™ 7 continues to improve workflow and make the job of Hearing Care Professionals faster and simpler than ever before. Combined with binax BestSound Technology, the latest-generation Connex 7 software helps you take full advantage of the new technology in just a few quick clicks. It supports the whole range of hearing aids – including support for the most advanced binaural aids on the market. Offer your customers optimal and individualized care by utilizing the new extended features, such as demonstrating binax' features or improved handling of microphone settings.

For more information, please consult our "Connex" brochure.



Understanding through demonstration

Your ability to explain changes to customers is a cornerstone of confidence in your services. Connex supports this by enabling use of the binaxDemo App. Now you can easily demonstrate changes you make in hearing aid settings to your customers.





Enhanced handling

Improved visualizations make it easier to see the effects of different microphone modes. Curve chart highlighting gives you direct visual feedback for gain adjustments you make. And the enlarged curve view lets you see all curves even more precisely.



More flexible tinnitus noiser function

In addition to the five static tinnitus therapy signals, the tinnitus noiser function now provides four ocean wave therapy signals, offering tinnitus patients a positive, soothing listening experience. Connex lets you easily configurate this flexible tinnitus therapy solution to meet any individual's needs.

| Products | Description | Item No. |
|---|---|----------|
| Connex | | |
|  | Connex 7, the powerful and easily learned fitting software for superior customer satisfaction, can be used to program the entire binax and micon range. | |
|  | Connex 7.4 DVD | 10675782 |
| | Connex 6.5 DVD | 10597368 |

Unity 3.



Enhance your workflow and
maximize customer satisfaction.



Siemens Unity™ 3 is a PC-controlled diagnostic and fitting system designed to simplify and speed up many of your daily work processes. Above and beyond that, it makes your business more cost-effective by minimizing downtime for calibration and service. The easy-to-use modular hardware and software components offer the utmost flexibility to fit your individual needs and its versatility enables faster and better measurements – so you can offer your customers a more pleasant, comfortable fitting experience with more precise results.



For more information, please consult our "Unity 3" brochure.

| Products | Description | Item No. |
|---|---|----------|
| Unity 3 Fitting Unit System | | |
|  | Unity 3 Fitting unit package incl. Unity 3 fitting unit, DC adapter, USB cord, client response, Unity quick installation guide, Unity software USB stick | 10823931 |
| Unity 3 Fitting Unit System | | |
|  | Unity 3 HIT package incl. coupler microphone, reference microphone, battery pills (S size), BTE tubes (3 pcs), ITE paddy, BTE coupler attachment, ITE coupler attachment, 2cc coupler attachment, body coupler attachment, USB cord, Hi-Pro extension cable, Unity quick installation guide, Unity software USB stick | 10823932 |

| Description | Item No. |
|------------------------------|----------|
| Transducers | |
| Real ear probe | 10823933 |
| EAR 3A | 10823934 |
| DD45 | 10823935 |
| Inserts earphone | 10823936 |
| HDA200 | 10823937 |
| Bone conductor – B71 | 10823938 |
| Accessories | |
| USB Isolation cord | 10823939 |
| Free field speaker (one) | 10824086 |
| Headset with mic (talk over) | 10823940 |
| Microphone (talk back) | 10823941 |
| Client response | 10823942 |
| Ear hanger (4 pcs) | 10823943 |
| DC adaptor | 10823944 |
| USB cord | 10824087 |
| Sound chamber cables | 10823946 |
| HIT ITE coupler attachment | 10823948 |
| HIT BTE coupler attachment | 10823949 |
| HIT 2cc coupler attachment | 10823947 |

| Description | Item No. |
|-------------------------------|----------|
| Accessories | |
| HIT body couple attachment | 10823950 |
| HIT coupler microphone | 10823951 |
| HIT reference microphone | 10823935 |
| HIT Hi-Pro extension cable | 10823936 |
| Consumables | |
| HIT ITE paddy | 08786084 |
| Insert tube with nipples | 10823956 |
| Probe tube guide (8 pcs) | 10823957 |
| Probe tubes (36 pcs) | 10823958 |
| BTE tube (3 pcs) | 10823959 |
| ITE foam | 10824088 |
| BTE foam | 10824089 |
| Battery Pills | |
| Battery pills (5 sizes) | 10823962 |
| Battery pill 5A (one piece) | 10823967 |
| Battery pill 10A (one piece) | 10823968 |
| Battery pill 312A (one piece) | 10823969 |
| Battery pill 13A (one piece) | 10823970 |
| Battery pill 675A (one piece) | 10823971 |

Workflow and Solutions.


| Products | Description | Item No. |
|--|---|----------------------------|
| Wireless Programming | | |
|  ConnexxLink | ConnexxLink wireless programming system for hearing aids equipped with e2e wireless capabilities | |
| | Package content: ConnexxLink wireless programming system, USB transmitter, power supply, lanyard, rechargeable battery, cradle | 10423664 |
| Programming | | |
|  Hi-Pro 2 USB NoahLink | Hi-Pro 2 USB PC Interface Hearing aid interface with Hi-Pro cable, using the serial USB port for communication between a PC and the Hi-Pro USB unit | 10297179 |
| | Hi-Pro 2 cable right | 02944986 |
| | Hi-Pro 2 cable left | 02944994 |
| | NoahLink incl. Bluetooth adapter NoahLink, wireless Bluetooth communication device to program hearing aids with ease, both in office and on location | 07784890 |
| | NoahLink cable | L: 07784908 R: 07784916 |
| Programming Adapter | | |
|  | Programming adapter, size 10A (1 pc) | 10824471 |
| | Programming adapter, size 312 (1 pc) | 10824469 |
| | Programming adapter, size 13 (1 pc) | 10824470 |
| | Programming adapter, Flex for custom (1 pc) | 03726994 |
| iScan | | |
|  | iScan | 08868007 |
| | <ul style="list-style-type: none"> Optimized scanning process: Up to twice as fast as the previous version (approx. 1.5 – 2 min) Open data format: The iScan data can also be used for placing orders with other manufacturers and otoplastic laboratories Hardware can be reused: Only software updates are necessary | |
| | Software | 10055228 |
| HearCheck | | |
|  | HearCheck™ Screener | 10175183 |
| | <ul style="list-style-type: none"> Check the hearing ability quickly and easily and determine whether the person has normal hearing or would require hearing aids Determine whether an open fitting would be suitable | |




eArena.

Acclimatization assistance for added value.

eArena™ is a useful auditory training resource for any business model. This 20-session interactive program helps patients increase listening competence when using their hearing aids in real life situations. Use it as a basis for acclimatization classes, give it to premium customers as an added value service, or share it with all your customers to increase customer loyalty and facilitate acceptance. Reduce the number of repeat visits and returns while increasing your opportunities to offer outstanding customer service with this proven tool.









| Products | Description | Item No. |
|---|--|----------|
| eArena | | |
|  | eArena auditory training for improved listening competence | 10826826 |

Point of Sale Support.

Along with a complete portfolio of highly competitive hearing instruments, we also supply Hearing Care Professionals with outstanding retail solutions. They guarantee an optimal brand and product presentation in the Hearing Care Professional's store. The range is based on the proven and popular product displays as well as new illuminated backlight advertising displays. These ensure a high-quality brand appearance at the point of sales. The comprehensive range of retail marketing materials provides everything consumers need for a satisfying purchasing experience.




Please contact your sales representative if you need detailed information about the retail marketing portfolio.



| Product Dummies | Description | Item | Item No. |
|---|--|--|--|
| Product Dummy Portfolio | | | |
|  | RIC binax dummies with magnet and dummy receiver. Tray packaging and dummy color sets also available. | 1. Ace binax dummy 2. Pure binax dummy 3. Carat binax dummy | A91SAT-02350-99G1-5000 A91SAT-02351-99G1-5000 A91SAT-02352-99G1-5000 |
|  | BTE binax dummies with magnet and dummy ThinTube (Motion SX) or dummy earhook (Motion PX). Tray packaging and dummy color sets also available. | 1. Motion SX binax dummy 2. Motion PX binax dummy | A91SAT-02463-99G1-5000 A91SAT-02465-99G1-5000 |
|  | ITE binax dummies with magnet. Tray packaging and ITE dummy set also available. | 1. Insio binax IIC dummy 2. Insio binax CIC dummy 3. Insio binax ITE dummy set, incl. CIC, ITC, ITC dir controls dummies | A91SAT-02470-99G1-5000 A91SAT-02468-99G1-5000 A91SAT-02372-99G1-5000 |
|  | Accessories – easyTek dummy with magnet. | easyTek dummy | A91SAT-02333-99G1-5000 |
| Dummies in Blister Pack | | | |
|  | Product dummy in blister pack. | 1. Ace binax 2. Pure binax 3. Carat binax 4. Motion SX binax 5. Motion PX binax | A91SAT-02350-99G1-5000 A91SAT-02351-99G1-5000 A91SAT-02352-99G1-5000 A91SAT-02473-99G1-5000 A91SAT-02474-99G1-5000 |
| Dummy Mailer | | | |
|  | RIC, BTE & ITE product dummies in blister packs. Can be used for mailings etc. | 1. RIC binax Dummy Mailer 2. Motion binax BTE Dummy Mailer 3. Insio binax ITE Dummy Mailer | A91SAT-02353-99G1-5000 A91SAT-02475-99G1-5000 A91SAT-02476-99G1-5000 |

| Product Presentation Displays | Description | Item | Item No. |
|---|---|---|--|
| Podium Displays | | | |
|  | Standalone podium display with alloy foot and metal bars in uprights. Suitable for all dummies with magnet. | 1. Small MPH for 2 dummies 2. Large MPH for 4 dummies 3. Portfolio MPH for 8 dummies | A91SAT-02299-99C1-5000 A91SAT-02300-99C1-5000 A91SAT-02301-99C1-5000 |
| Product Holders | | | |
|  | Single Product Holder with special foot including illumination. | 1. Illuminated Single Product Holder | A91SAT-01887-99C1-5000 |
| | Single Product Holder with foot. | 2. Single Product Holder, orange 3. Single Product Holder, aqua | A91SAT-01888-99C1-5000 A91SAT-01889-99C1-5000 |
| Product Display Systems | | | |
|  | Large display including 15" display and 8 Single Product Holders (H 502.6 mm, W 427 mm, D 424 mm). | 1. Digital Large Display | A91SAT-00905-99C2-5000 |
| | Small display including 3 Single Product Holders (H 285.2 mm, W 212 mm, D 246.6 mm). | 2. Small Display | A91SAT-00909-99C2-5000 |
| Counseling Tools | | | |
| Consultation Box | | | |
|  | Including 3 empty product trays, 1 accessory drawer, 1 real ear plate, 1 removable lid tool with mirror (H 143 mm, W 242 mm, D 195 mm, closed). | Consultation Box | A91SAT-00909-99C2-5000 |
| Presentation Trays | | | |
|  | RIC, BTE, ITE presentation trays including product dummies and color plates or print insert showing available products & colors. Suitable for Consultation Box. | 1. RIC binax Presentation tray 2. Motion binax BTE Presentation tray 3. Insio binax ITE Presentation tray | A91SAT-02348-99G1-5000 A91SAT-02477-99G1-5000 A91SAT-02478-99G1-5000 |





Point of Sale Support.





| In-Store Advertising Displays | Description | Item | Item No. |
|--|--|--|------------------------|
| Billboard Frames | | | |
| Billboard Floorstanding Display | | | |
|  | Billboard Floorstanding Displays fitting poster snapframe or lightbox in DIN A3 or banner format and leaflet holders. Double sided configuration on request. | 1. Billboard Floorstanding Banner Lightbox | A91SAT-02188-99C1-5000 |
| | | 2. Billboard Floorstanding Banner Snapframe | A91SAT-02189-99C1-5000 |
| | | 3. Billboard Floorstanding A3 Lightbox, 1 Leaflet Holder | A91SAT-02191-99C1-5000 |
| | | 4. Billboard Floorstanding A3 Snapframe, 2 Leaflet Holders | A91SAT-02320-99C1-5000 |
| Billboard Counter Display | | | |
|  | Billboard Counter Displays fitting poster snapframe or lightbox in DIN A3 or used as leaflet holder | 1. Billboard Counter A3 Lightbox | A91SAT-02192-99C1-5000 |
| | | 2. Billboard Counter A3 Snapframe | A91SAT-02193-99C1-5000 |
| | | 3. Billboard Counter Leaflet Holder | A91SAT-02194-99C1-5000 |
| Billboard Standalone Frames | | | |
|  | Billboard poster snapframe or lightbox in DIN A3 or banner format with wall-mounting kit | 1. Billboard Banner Lightbox | A91SAT-02195-99C1-5000 |
| | | 2. Billboard Banner Snapframe | A91SAT-02263-99C1-5000 |
| | | 3. Billboard A3 Lightbox | A91SAT-02264-99C1-5000 |
| | | 4. Billboard A3 Snapframe | A91SAT-02265-99C1-5000 |

Jewel Cases.

Add an element of luxury to the customer experience by providing attractive cases for purchased hearing aids. The basic and value options are sturdy packages that highlight the well-crafted technology inside. And the premium options add a sense of elegance.



| In-Store Advertising | Description | Item | Item No. |
|---|--|----------------------------|------------------------|
| Driver Pack – Promotion & Decoration Package for point of sales | | | |
| Driver Packages “Enjoy the sound of Seasons” | | | |
|  | | | |
| Driver Pack Winter Edition “Enjoy the sound of Winter” | | | |
|  | Seasonal decoration packages for in-store advertising. Easy-to-adapt POS materials. Seasonal edition available in digital format for local production. | 1. Banner | A91SAT-02437-99C1-7600 |
| | | 2. Poster | A91SAT-02438-99C1-7600 |
| | | 3. Convex Frame Inserts | A91SAT-02439-99C1-7600 |
| | | 4. Counter Display | A91SAT-02440-99C1-7600 |
| | | 5. Window Sticker | A91SAT-02441-99C1-7600 |
| | | 6. Cube big | A91SAT-02442-99C1-7600 |
| | | 7. Cube small | A91SAT-02443-99C1-7600 |
| Standard Frames | | | |
| Poster Clip Frame | | | |
|  | Poster Clip Frame for DIN A1 portrait posters, with Siemens branding and various hanging possibilities | Poster Clip Frame | A91SAT-01730-99C1-5000 |
| Convex Frames | | | |
|  | Convex Frames can hold two print inserts on either side and are easily updatable. Specific print inserts available for all products. Set of 2 prints | Convex Frame Poster Holder | A91SAT-00667-99C1-7600 |






| Products | Description | Item No. | Products | Description | Item No. |
|---|-------------|----------|---|-------------|----------|
| Premium Jewel Case | | | Premium Box | | |
|  | Large case | 10672623 |  | One size | 10672625 |
| | Small case | 10672624 | | | |
| Basic Case | | | Value Box | | |
|  | One size | 10821648 |  | One size | 10672626 |

All Sizes at a Glance – RIC and BTE.

Actual sizes in cm



Detailed Feature Overview RIC.

| Type | Ace binax | Pure binax | Carat/Carat A binax | Pure micon | Orion RIC |
|---|---|---|---|---|---|
| Performance levels | 7bx/5bx/3bx New! | | | 7mi/5mi/3mi | Essential |
| Picture |  |  |  |  |  |
| Product details | | | | | |
| Matrix (2 ccm) | RIC S 45/108 dB M 55/113 dB P 60/118 dB HP 65/123 dB | S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB | S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB | S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB | S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB |
| Battery size | 10A | 312 | 13 | 312 | 312 |
| Power consumption | 0.9 – 1.0 mA | 0.9 – 1.1 mA | 0.9 – 1.1 mA | 0.8 – 1.1 mA | 0.8 – 1.1 mA |
| Features | | | | | |
| BestSound Technology | binax | binax | binax | micon | - |
| Frequency channels/ Compression controls | 7bx – 48/20 5bx – 32/16 3bx – 24/12 - | 7bx – 48/20 5bx – 32/16 3bx – 24/12 - | 7bx – 48/20 5bx – 32/16 3bx – 24/12 - | 7mi – 48/20 5mi – 32/16 3mi – 24/12 - | 16/8 |
| Hearing programs | 7bx – 6 5bx – 6 3bx – 6 - | 7bx – 6 5bx – 6 3bx – 6 - | 7bx – 6 5bx – 6 3bx – 6 - | 7mi – 6 5mi – 6 3mi – 6 - | 4 |
| Rechargeable | - | ● ¹⁶ | ● ^{4,22} | ● ¹⁶ | - |
| Volume control | - | ○ ¹ | ● | ○ ¹ | ○ ¹ |
| Push button | ● ¹ | ○ ¹ | - | ○ ¹ | ○ ¹ |
| T-coil | - | ● | ● | ● | ● |
| Battery door on/off function | ● | ● | ● | ● | ● |
| Audio input | - | ○ ³ | ○ ^{3,5} | ○ ⁶ | - |
| Tinnitus noiser | ● | ● | ● | ● | ● |
| Modifications | | | | | |
| Color conversion kit | ○ | ● | ● | ● | ● |
| Accessories | | | | | |
| Audio streamers | - | ● | ● | ● | - |
| Remote controls | - | ● | ● | ● | ● |
| Programming adapter/cable | size 10 | size 312 | size 13 | size 312 | size 312 |
| Fitting set | miniReceiver 2.0 | miniReceiver 2.0 | miniReceiver 2.0 | miniReceiver 2.0 | miniReceiver 2.0 |
| ConnexLink | - | ● | ● | ● | - |

(1) Exchangeable controls – programmable via Connexx
(3) With easyTek
(4) Carat binax models only

(5) Carat A binax models only
(6) With miniTek
(16) Only recommended with miniReceiver S
(22) Not recommended with miniReceiver HP









● Standard
○ Optional







Detailed Feature Overview BTE.

| Type | Motion SX binax | Motion SA binax | Motion PX binax | Siemens Life micon | Motion SX micon | Motion M micon | |
|---|---|---|---|--|---|---|------------------------|
| Performance levels | 7bx/5bx/3bx | | | 7mi/5mi/3mi | 7mi/5mi/3mi | | |
| Picture |  |  |  |  |  |  | |
| Product details | | | | | | | |
| Matrix (2ccm) | ThinTube/LifeTube Earhook | 53/125 dB 60/130 dB | 53/125 dB 60/130 dB | 65/127 dB 75/135 dB | 47/121 dB 57/123 dB | 53/125 dB 60/130 dB | 60/124 dB 65/130 dB |
| Battery size | 13 | 13 | 13 | 312 | 13 | 13 | |
| Power consumption | 1.0 mA | 1.0 mA | 1.1 – 1.3 mA | 0.8 mA | 1.0 mA | 1.1 mA | |
| Features | | | | | | | |
| BestSound Technology | binax | binax | binax | micon | micon | micon | |
| Frequency channels/ Compression controls | 7bx – 48/20 5bx – 32/16 3bx – 24/12 - | 7bx – 48/20 5bx – 32/16 3bx – 24/12 - | 7bx – 48/20 5bx – 32/16 3bx – 24/12 - | 7mi – 48/20 5mi – 32/16 3mi – 24/12 - | 7mi – 48/20 5mi – 32/16 3mi – 24/12 - | 7mi – 48/20 5mi – 32/16 3mi – 24/12 - | |
| Hearing programs | 7bx – 6 5bx – 6 3bx – 6 - | 7bx – 6 5bx – 6 3bx – 6 - | 7bx – 6 5bx – 6 3bx – 6 - | 7mi – 6 5mi – 6 3mi – 6 - | 7mi – 6 5mi – 6 3mi – 6 - | 7mi – 6 5mi – 6 3mi – 6 - | |
| Rechargeable | ● | - | ● | - | ● | - | |
| Volume control | ● ⁷ | ● ⁷ | ● | ○ ⁷ | ○ ⁷ | ● | |
| Push button | - | ● | ● | ● | ○ ⁹ | ● | |
| T-coil | ● | ● | ● | ○ ⁶ | ● | ● | |
| Battery door on/off function | ● | ● | - | ● | ● | ● | |
| LifeTube/Standard earhook | ● | ● | ● | ● | ● | ● | |
| Audio input | ○ ³ | ○ ⁸ | ○ ³ | ○ ⁶ | ○ ⁶ | ○ ⁸ | |
| Tinnitus noiser | ● | ● | ● | ● | ● | ● | |
| Modifications | | | | | | | |
| Color conversion kit | ● | ● | ● | ● | ● | ● | |
| Eyeglass adapter | - | - | ○ | - | ○ | ○ | |
| Accessories | | | | | | | |
| Audio streamers | ● | ● | ● | ● | ● | ● | |
| Remote controls | ● | ● | ● | ● | ● | ● | |
| Programming adapter/cable | size 13 | size 13 | size 13 | size 312 | size 13 | size 13 | |
| Fitting set | ThinTube | ThinTube | S-LifeTube | LifeTube | LifeTube | S-LifeTube | |
| ConnexLink | ● | ● | ● | ● | ● | ● | |

- (1) Exchangeable controls – programmable via Connexx
 - (2) Programmable via Connexx
 - (3) With easyTek
 - (4) With battery compartment modification
 - (5) With remote control
 - (6) With miniTek
 - (7) Programmable via Connexx
 - (8) With battery compartment modification
 - (9) With remote control
 - (10) With earhook and eyeglass modification
- Standard
○ Optional






| Motion P micon | Motion PX micon | Aquaris micon | Nitro micon | Orion/Sirion S | Orion/Sirion M | Orion/Sirion P | Orion SP |
|---|---|---|---|---|---|---|---|
| 7mi/5mi/3mi | | | 7mi/5mi/3mi | 7mi/3mi | Essential | | |
|  |  |  |  |  |  |  |  |
| 65/127 dB 75/135 dB | 65/127 dB 75/135 dB | 49/122 dB 65/129 dB | - 82/141 dB | 45/124 dB 55/124 dB | 53/125 dB 60/130 dB | 65/125 dB 70/130 dB | - 82/140 dB |
| 13 | 13 | 13 | 675 | 312 | 13 | 13 | 675 |
| 1.1 – 1.3 mA | 1.1 – 1.3 mA | 1.0 – 1.1 mA | 2.1 mA | 0.8 mA | 1.0 mA | 1.0 – 1.3 mA | 2.1 mA |
| micon | micon | micon | micon | - | - | - | - |
| 7mi – 48/20 5mi – 32/16 3mi – 24/12 - | 7mi – 48/20 5mi – 32/16 3mi – 24/12 - | 7mi – 48/20 5mi – 32/16 3mi – 24/12 - | 7mi – 24/16 - 3mi – 24/12 - | 16/8 12/6 | 16/8 12/6 | 16/8 12/6 | 16/8 |
| 7mi – 6 5mi – 6 3mi – 6 - | 7mi – 6 5mi – 6 3mi – 6 - | 7mi – 6 5mi – 6 3mi – 6 - | 7mi – 6 - 3mi – 6 - | 4 | 4 | 4 | 4 |
| - | ● | - | - | - | - | - | - |
| ● | ● | ○ ⁹ | ● | ○ ⁷ | ○ ⁷ | ● | ● |
| ● | ● | ○ ⁹ | ● ⁷ | ● ⁷ | ○ ⁷ | ● | ● |
| ● | ● | ○ ⁶ | ● ⁷ | - | ● | ● | ● |
| ● | ● | ● | ● | ● | ● | ● | ● |
| ● | ● | ● | ● | ● | ● | ● | ● |
| ○ ⁸ | ○ ⁶ | ○ ⁶ | ○ ⁸ | - | ○ ⁸ | ○ ⁸ | ○ ⁸ |
| ● | ● | ● | ● | ● | - | ● | ● |
| ● | ● | ● | only top housing | ● | ● | ● | - |
| ○ | ○ | - | - | - | - | - | - |
| ● | ● | ● | ● | - | - | - | - |
| ● | ● | ● | ● | ● | - | ● | ● |
| size 13 | size 13 | CS 44 | CS 44 | size 312 | size 312 | size 312 | CS 44 |
| S-LifeTube | S-LifeTube | LifeTube | - | LifeTube | LifeTube | LifeTube | - |
| ● | ● | ● | ● | ● | ● | ● | ● |

Detailed Feature Overview BTE.

| Type | Intuis Life | Intuis S Dir | | Intuis Dir | | Intuis SP Dir | | |
|---|---|------------------------|---|----------------|--|----------------|---|--|
| Performance levels | | | | | | | | |
| Basic | | | | | | | | |
| Picture |  | |  | |  | |  | |
| Product details | | Pro S Dir | S Dir | Pro Dir | Dir | Pro SP | Dir | |
| Matrix (2 ccm) | LifeTube Earhook | 35/115 dB 45/117 dB | 53/123 dB 55/124 dB | - 60/130 dB | | - 80/138 dB | | |
| Battery size | | 312 | 13 | 13 | | 675 | | |
| Power consumption | | 0.5 mA | 0.7 mA | 0.8 mA | | 1.0 mA | | |
| Features | | | | | | | | |
| Frequency channels/ Compression controls | | 4/4 | 4/4 | 4/4 | | 4/4 | | |
| Trimmers | | - | - | - | | - | | |
| Hearing programs | | 1 | 1-4 4 | 1-4 4 | 4 | 1-4 4 | 4 | |
| Rechargeable | | - | - | - | | - | | |
| Volume control | | - | - | ● | | ● | | |
| Push button | | - | ● | ● | | ● | | |
| T-coil | | - | ● | ● | | ● | | |
| Battery door on/off function | | ● | ● | ● | | ● | | |
| LifeTube/Standard earhook | | ● | ● | ● | | ● | | |
| Audio input | | - | - | ● | | ● | | |
| Tinnitus noiser | | - | - | - | | - | | |
| Modifications | | | | | | | | |
| CROS | | ● ¹³ | ● ¹³ | ● | | ● | | |
| BICROS | | ● ¹³ | ● ¹³ | ● | | ● | | |
| Eyeglass adapter | | ○ | ○ | ○ | | ○ | | |
| Accessories | | | | | | | | |
| Programming adapter/cable | | CS 44 | CS 44 | CS 44 | | CS 44 | | |
| Fitting set | | LifeTube | LifeTube | - | | - | | |

(13) With eyeglass modification

● Standard
○ Optional

| Lotus Pro 2 M | Lotus Pro 2 SP | Lotus M | | Lotus P | | Lotus SP | | | |
|---|----------------|---|----------------|---|----------------|---|-------------|---|--|
| Basic | | | | | | | | | |
| Basic | | | | Basic | | | | | |
|  | |  | |  | |  | |  | |
| 2M | 2SP | Pro | 23 | Pro | 23 | Pro | 23 | | |
| - 60/130 dB | - 80/138 dB | 47/126 dB 60/130 dB | - 60/130 dB | - 70/133 dB | - 70/132 dB | - 80/138 dB | | | |
| 13 | 675 | 13 | | 13 | | 675 | | | |
| 0.8 mA | 1.6 mA | 0.8-1.0 mA | 1.0 mA | 0.7 mA | 0.9 mA | 2.0 mA | 0.9 mA | | |
| Features | | | | | | | | | |
| 2/1 | 2/1 | 2/1 | 2 | 2/1 | 2 | 2/1 | 2 | | |
| - | - | - | NH/MPO/Gain | - | NH/MPO/Gain | - | NH/MPO/Gain | | |
| 3+A | 3+A | 3+A | 2 | 2+A | 2 | 2+A | 2 | | |
| - | - | - | - | - | - | - | - | | |
| ● | ● | ● | ● | ● | ● | ● | ● | | |
| ● | ● | ● | ● | ● | ● | ● | ● | | |
| ● | ● | ● | ● | ● | ● | ● | ● | | |
| ● | ● | ● | ● | ● | ● | ● | ● | | |
| ● | ● | ● | ● | ● | ● | ● | ● | | |
| ● | ● | ● | ● | ● | ● | ● | ● | | |
| ● | ● | ○ | - | ● | - | ● | - | | |
| Accessories | | | | | | | | | |
| CS 44 | CS 44 | CS 44 | | CS 44 | | CS 44 | | | |
| - | - | LifeTube | - | - | - | - | | | |

Detailed Feature Overview Custom.


| Type | Insio binax IIC / CIC | Insio binax ITC | Insio binax ITE | Insio micon IIC/CIC | Insio micon ITC | Insio micon ITE |
|---|---|---|---|---|---|---|
| Performance levels | 7bx/5bx/3bx | | | 7mi/5mi/3mi | | |
| Picture | | | | | | |
| Product details | | | | | | |
| Matrix (2 ccm) | 50/113 ¹¹ dB | 55/118 ¹¹ dB | 60/123 ¹¹ dB | 50/113 ¹¹ dB | 55/118 ¹¹ dB | 60/123 ¹¹ dB |
| Battery size | 10 | 312/10 | 13/312 | 10 | 312/10 | 13/312 |
| Power consumption | 0.9 – 1.0 mA | 1.1 mA | 1.1 mA | 0.9 mA | 1.0 mA | 1.0 mA |
| Features | | | | | | |
| BestSound Technology | binax | binax | binax | micon | micon | micon |
| Frequency channels/ Compression controls | 7bx – 48/20 5bx – 32/16 3bx – 24/12 | 7bx – 48/20 5bx – 32/16 3bx – 24/12 | 7bx – 48/20 5bx – 32/16 3bx – 24/12 | 7mi – 48/20 5mi – 32/16 3mi – 24/12 | 7mi – 48/20 5mi – 32/16 3mi – 24/12 | 7mi – 48/20 5mi – 32/16 3mi – 24/12 |
| Hearing programs | 7bx – 6 5bx – 6 3bx – 6 | 7bx – 6 5bx – 6 3bx – 6 | 7bx – 6 5bx – 6 3bx – 6 | 7mi – 6 ¹² 5mi – 6 ¹² 3mi – 6 ¹² | 7mi – 6 5mi – 6 3mi – 6 | 7mi – 6 5mi – 6 3mi – 6 |
| Volume control | ○ ^{6,12} | ○ | ○ | ○ ^{6,12} | ○ | ○ |
| Push button | ○ ¹² | ○ | ○ | ○ ¹² | ○ | ○ |
| T-coil | - | ○ | ○ | - | ○ | ○ |
| Audio input | ○ ^{3,12} | ○ ³ | ○ ³ | ○ ^{6,12} | ○ ⁶ | ○ ⁶ |
| Tinnitus noiser | ● ¹² | ● | ● | ● ¹² | ● | ● |
| Optivent | ● | ● | ● | ○ | ● | ● |
| Accessories | | | | | | |
| Audio streamers | ● ¹² | ● | ● | ● ¹² | ● | ● |
| Remote controls | ● ¹² | ● | ● | ● ¹² | ● | ● |
| Programming adapter | Flex connector | Flex connector | Flex connector | Flex connector | Flex connector | Flex connector |

(3) With easyTek (6) With miniTek (11) Max. values (12) Not available for IIC


● Standard ○ Optional

| Nitro micon CIC | Nitro micon ITC | Nitro micon ITE | Orion/Sirion CIC | | Orion/Sirion ITC | | Orion/Sirion ITE | |
|---|---|---|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|-------------------------|
| 7mi/3mi | | | Essential | | | | | |
| | | | | | | | | |
| for all models from 55/118 ¹¹ dB up to 71/130 ¹¹ dB | for all models from 55/118 ¹¹ dB up to 71/130 ¹¹ dB | for all models from 55/118 ¹¹ dB up to 71/130 ¹¹ dB | Orion | Sirion | Orion | Sirion | Orion | Sirion |
| 50/113 ¹¹ dB | 50/113 ¹¹ dB | 50/113 ¹¹ dB | 55/118 ¹¹ dB | 55/118 ¹¹ dB | 55/118 ¹¹ dB | 55/118 ¹¹ dB | 60/123 ¹¹ dB | 60/123 ¹¹ dB |
| 10 | 312/10 | 13/312 | 10 | 312/10 | 10 | 312/10 | 13/312 | 13/312 |
| 1.0 mA | 1.0 – 1.2 mA | 1.0 – 1.2 mA | 0.9 mA | 1.0 mA | 0.9 mA | 1.0 mA | 1.0 mA | 1.0 mA |
| micon | micon | micon | - | - | - | - | - | - |
| 7mi – 48/20 | 7mi – 48/20 | 7mi – 48/20 | 16/8 | 12/6 | 16/8 | 12/6 | 16/8 | 12/6 |
| - | - | - | - | - | - | - | - | - |
| 3mi – 24/12 | 3mi – 24/12 | 3mi – 24/12 | - | - | - | - | - | - |
| - | - | - | - | - | - | - | - | - |
| 7mi – 6 | 7mi – 6 | 7mi – 6 | 4 | 4 | 4 | 4 | 4 | 4 |
| - | - | - | - | - | - | - | - | - |
| 3mi – 6 | 3mi – 6 | 3mi – 6 | - | - | - | - | - | - |
| - | - | - | - | - | - | - | - | - |
| ○ | ○ | ○ | - | ○ | - | ○ | - | ○ |
| ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ | ○ |
| - | ○ | ○ | - | ○ | - | ○ | - | ○ |
| ○ ⁶ | ○ ⁶ | ○ ⁶ | - | - | - | - | - | - |
| ● | ● | ● | ● | - | ● | - | ● | - |
| ● | ● | ● | ● | ● | ● | ● | ● | ● |
| - | - | - | - | - | - | - | - | - |
| - | - | - | ● | - | ● | - | ● | - |
| Flex connector | Flex connector | Flex connector | Flex connector | Flex connector | Flex connector | Flex connector | Flex connector | Flex connector |

Custom instruments features



Custom instruments









Standard skin colors for custom instruments

| Description | Type | Insio | Nitro | Orion/Sirion |
|------------------|---------------|-------|-------|--------------|
| Beige | ITC, ITE | ● | ● | ● |
| Tan | ITC, ITE | ● | ● | ● |
| Mocha | CIC, ITC, ITE | ● | ● | ● |
| Brown | CIC, ITC, ITE | ● | ● | ● |
| Dark Brown | ITC, ITE | ● | ● | ● |
| White/Dark Brown | IIC | ● | | |






Detailed Feature Overview Custom.

| Type | Intuis CIC | Intuis ITC | Intuis ITE | Lotus CIC | Lotus ITC | Lotus ITE | |
|---|---|---|---|---|---|--|-------------------------|
| Performance levels | Basic | | | Basic | | | |
| Picture |  |  |  |  |  |  | |
| Product details | | | | Pro | 23 | Pro | 23 |
| Matrix (2 ccm) | 47/113 ¹¹ dB | 50/118 ¹¹ dB | 60/123 ¹¹ dB | 47/113 ¹¹ dB | 40/113 ¹¹ dB | 50/118 ¹¹ dB | 60/123 ¹¹ dB |
| Battery size | 10 | 312 | 13 | 10 | 312 | 13 | |
| Power consumption | 0.6 mA | 0.6 mA | 0.6 mA | 0.6 mA | 0.7 mA | 0.6 mA | 0.5 mA |
| Features | | | | | | | |
| Frequency channels/ Compression controls | 4/4 | 4/4 | 4/4 | 2/1 | 2 | 2/1 | 2 |
| Trimmers | - | - | - | - | NH/MPO/NR | - | NH/MPO/NR |
| Hearing programs | ○ | ○ | ○ | ○ | 1 | ○ | 1 |
| Volume control | - | ○ | ○ | - | Manual VC ¹⁴ | ○ | - |
| Push button | ○ | ○ | ○ | - | - | ○ | - |
| T-coil | - | ○ | ○ | - | - | ○ | ○ |
| Audio input | - | - | - | - | - | - | - |
| Tinnitus noiser | - | - | - | - | - | - | - |
| Accessories | | | | | | | |
| Programming adapter | Flex connector | Flex connector | Flex connector | Flex connector | Flex connector | Flex connector | Flex connector |

(11) Max. values
 (14) With screw set volume control for CIC

● Standard
 ○ Optional

Standard Skin Colors for Custom Instruments

| Description | Type | Intuis | Lotus |
|--|---------------|--------|-------|
|  Beige | ITC, ITE | ● | ● |
|  Tan | ITC, ITE | ● | ● |
|  Brown | CIC, ITC, ITE | ● | ● |



Color Overview RIC and BTE.

| Description | Color No. | Ace binax | Pure binax | Carat binax | Motion binax | Pure micon | Siemens Life micon |
|----------------|-----------|-----------|------------|-------------|-----------------|------------|--------------------|
| Beige | 1 | ● | ● | ● | ● | ● | ● |
| Granite | 2 | ● | ● | ● | ● | ● | ● |
| Grey | 3 | ● | ● | ● | ● | ● | ● |
| Brown | 4 | | ● | ● | ● | ● | ● |
| Black | 5 | ● | ● | ● | ● | ● | ● |
| Silver | 6 | ● | ● | ● | ● | ● | ● |
| Pearl White | 7 | ● | ● | ● | ● ²¹ | ● | ● |
| Golden Blonde | 30 | ● | ● | ● | ● | ● | ● |
| Sandy Brown | 31 | ● | ● | ● | ● | ● | ● |
| Dark Champagne | 37 | | ● | ● | ● | ● | ● |
| Sandy | 51 | | | | | | |
| Silver Lake | 53 | | | | | | |
| Seagull | 54 | | | | | | |
| Stony Beach | 55 | | | | | | |
| Aqua | 56 | | | | | | |
| Candy Pink | 60 | | | | ● ¹⁹ | | |
| Galactic Blue | 61 | | | | ● ¹⁹ | | |
| Dark Granite | 66 | ● | ● | ● | ● | ● | ● |
| Spirit | C5 | ● | ● | ● | ● ²⁰ | ● | |
| Sparkling | C10 | | | | | | ● |
| Elegance | C12 | | | | ● ²⁰ | | |

(17) Only Motion M and P micon
 (18) Only Motion SX and PX micon
 (19) Only Motion SA binax
 (20) Only Motion SA and SX binax
 (21) Only Motion PX binax

| Motion micon | Aquaris micon | Nitro micon | Orion RIC | Orion/Sirion S/M/P | Orion SP | Intuis | Intuis Pro | Lotus Pro | Lotus Trimmer |
|-----------------|---------------|-------------|-----------|--------------------|----------|--------|------------|-----------|---------------|
| ● | | ● | ● | ● | ● | ● | ● | ● | ● |
| ● | | ● | ● | | ● | ● | ● | | |
| ● | | ● | ● | ● | ● | ● | ● | ● | ● |
| ● ¹⁸ | | ● | ● | ● | ● | ● | ● | ● | ● |
| ● | | ● | ● | ● | ● | ● | | | |
| ● | | | ● | ● | | ● | | | |
| ● ¹⁸ | | | | | | | | | |
| ● ¹⁸ | | | | | | | | | |
| ● | | | | | | | | | |
| ● ¹⁸ | | | | | | | | | |
| | ● | | | | | | | | |
| | ● | | | | | | | | |
| | ● | | | | | | | | |
| | ● | | | | | | | | |
| | ● | | | | | | | | |
| ● ¹⁷ | | | | | | | | | |
| ● ¹⁷ | | | | | | | | | |
| ● | | ● | ● | ● | ● | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |



Color Flipper, version October 2014
 Presenting RIC, BTE and ITE colors, and product-specific overview of available colors per product.

Item No. A91SAT-02349-99C1-7600

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

* Two clinical studies have shown that binax provides better than normal hearing in certain demanding environments (University of Northern Colorado, 2014; Oldenburg Hörzentrum, 2013): Speech Reception Thresholds (SRT) in cocktail-party situations improved up to 2.9 dB for wearers with mild to moderate hearing loss using Carat binax or Pure binax hearing aids with narrow directionality, compared to people with normal hearing.

Android and Google Play are trademarks of Google Inc. Apple App Store is a trademark of Apple Inc.

The Bluetooth word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos GmbH is under license. Other trademarks and trade names are those of their respective owners.

GORE and designs are trademarks of W. L. Gore & Associates.

Legal Manufacturer
Sivantos GmbH
Henri-Dunant-Strasse 100
91058 Erlangen
Germany

Sivantos GmbH is a Trademark Licensee of Siemens AG.