

NEW

PerfectClean® works with a refill

Antibacterial solution disinfects and removes ear wax.
Effect stressed with the Vortex Power.



PerfectClean®
given to the user



PerfectClean®

The first system to clean,
dry and disinfect hearing aids



Customer Tests show

Customers tests (B.Chaix) show that after use,
95% of customers consider PerfectClean® to be very efficient
and easy to use.

70% of customers would recommend the PerfectClean® to
their friends

100% of users agree about the High-End Design!

Produced by :



MG Development America Inc.
187 36th Street
Pittsburgh, PA 15201 – USA
www.mg-dev.com

Mg Development is dedicated to the development and the
production of **Medical Devices** for the maintenance of Hearing
Systems and the Ear Care.



Hear well every day

PerfectClean®

The first system to clean, disinfect, and dry hearing aids with a press of a button

What are the BENEFITS of the PerfectClean®?

- Effectively removes earwax from hearing aid devices
- Helps to prevent irritation
- Maintains optimal performance of hearing aid devices
- Easy to Use

1

Place hearing aids inside

Fits all types (standard and custom)

2


Press button

3 actions in 1h30

Cleaning: earwax removal solution

Disinfecting: active anti-bacterial agent + UV-C system

Drying: forced-air fan drying system



Vortex cleaning action eliminates earwax



Forced-air fan drying system

- Temperature self-regulated at 45°C +/-3°C
- Controlled ventilation
- Optimized and consistent air flow
- Extraction of moisture



UV-C light completely sanitizes hearing aids



- Antibacterial solution disinfects and removes ear wax
- Easy to refill
- Visual indication of the solution level

PerfectClean®